



**Physical Characteristics:**

particle size: 5µm, spherical  
pore size: 60Å  
carbon load: 27%  
endcap: fully endcapped  
pH range: 2.5 to 8  
temperature limit: 80°C



.....	Page #
herbicides .....	.712
pesticides .....	.716

**Allure® Columns: 5µm particles; 60Å pore size**

Small pore size in a high-purity, Type B silica provides a large surface area. High carbon loads, highly retentive. An excellent choice for evaporative light scattering (ELSD) and MS detectors, in which more organic solvent in the mobile phase gives better sensitivity.

**Allure® C18 Columns (USP L1)**

**Chromatographic Properties:**

Most retentive of our alkyl stationary phases due to large surface area of the base silica and high-density bondings. Provides excellent peak shapes for a wide range of compounds.

Length	1.0mm ID		2.1mm ID		3.2mm ID		4.6mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
<b>5µm Columns</b>								
30mm	9164531		9164532		9164533		9164535	
50mm	9164551		9164552		9164553		9164555	
100mm	9164511		9164512		9164513		9164515	
150mm	9164561		9164562		9164563		9164565	
200mm	9164521		9164522		9164523		9164525	
250mm	9164571		9164572		9164573		9164575	

**Allure® C18 Guard Cartridges**

Guard Cartridges	3-pk. (10 x 2.1mm)	3-pk. (10 x 4.0mm)	2-pk. (20 x 2.1mm)	2-pk. (20 x 4.0mm)	price
Allure C18 Guard Cartridge	916450212	916450210	916450222	916450220	\$148

**Allure® C18 HPLC Prep Columns**

Length	10mm ID		21.2mm ID		30mm ID		50mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
<b>5µm Columns</b>								
50mm	9164557		9164558		9164559		9164550	
100mm	9164517		9164518		9164519		9164510	
150mm	9164567		9164568		9164569		9164560	
250mm	9164577		9164578		9164579		9164570	

**Physical Characteristics:**

particle size: 5µm, spherical  
pore size: 60Å  
carbon load: 12%  
endcap: fully endcapped  
pH range: 2.5 to 8  
temperature limit: 80°C



.....	Page #
herbicides .....	.714

**Allure® Basix Columns (USP L10)**

**Chromatographic Properties:**

Highly retentive propyl cyano phase. Excellent choice for basic compounds and for analytes containing amine group functionality.

Length	1.0mm ID		2.1mm ID		3.2mm ID		4.6mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
<b>5µm Columns</b>								
30mm	9161531		9161532		9161533		9161535	
50mm	9161551		9161552		9161553		9161555	
100mm	9161511		9161512		9161513		9161515	
150mm	9161561		9161562		9161563		9161565	
200mm	9161521		9161522		9161523		9161525	
250mm	9161571		9161572		9161573		9161575	

**Allure® Basix Guard Cartridges**

Guard Cartridges	3-pk. (10 x 2.1mm)	3-pk. (10 x 4.0mm)	2-pk. (20 x 2.1mm)	2-pk. (20 x 4.0mm)	price
Allure Basix Guard Cartridge	916150212	916150210	916150222	916150220	

**Allure® Basix HPLC Prep Columns**

Length	10mm ID		21.2mm ID		30mm ID		50mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
<b>5µm Columns</b>								
50mm	9161557		9161558		9161559		9161550	
100mm	9161517		9161518		9161519		9161510	
150mm	9161567		9161568		9161569		9161560	
250mm	9161577		9161578		9161579		9161570	



also  
**available**

**Bulk Packing Materials**  
See page 194.