

Ultra PFP Columns (USP L43)

Chromatographic Properties:

A pentafluorophenyl phase. Unique selectivity for compounds containing organohalogens or other basic functional groups.

Length	1.0mm ID		2.1mm ID		3.2mm ID		4.6mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
3µm Columns								
30mm	9176331		9176332		9176333		9176335	
50mm	9176351		9176352		9176353		9176355	
100mm	9176311		9176312		9176313		9176315	
5µm Columns								
30mm	9176531		9176532		9176533		9176535	
50mm	9176551		9176552		9176553		9176555	
100mm	9176511		9176512		9176513		9176515	
150mm	9176561		9176562		9176563		9176565	
200mm	9176521		9176522		9176523		9176525	
250mm	9176571		9176572		9176573		9176575	

3.0mm ID available on request for 3µm particle applications above 4,000 psi (275 Bar).

Ultra PFP Guard Cartridges

Guard Cartridges	3-pk. (10 x 2.1mm)	3-pk. (10 x 4.0mm)	2-pk. (20 x 2.1mm)	2-pk. (20 x 4.0mm)	price
Ultra PFP Guard Cartridge	917650212	917650210	917650222	917650220	

Ultra PFP HPLC Prep Columns

Length	10mm ID		21.2mm ID		30mm ID		50mm ID	
	cat.#	price	cat.#	price	cat.#	price	cat.#	price
5µm Columns								
50mm	9176557		9176558		9176559		9176550	
100mm	9176517		9176518		9176519		9176510	
150mm	9176567		9176568		9176569		9176560	
250mm	9176577		9176578		9176579		9176570	

Physical Characteristics:

particle size: 3µm or 5µm, spherical
pore size: 100Å
carbon load: 7%
endcap: fully endcapped
pH range: 2.5 to 8
temperature limit: 80°C



.....	Page #
purines, pyrimidines	745




ordering note

To order a 2.1mm, 3.2mm, or 4.6mm ID column with a Trident Integral Inlet Fitting, add "-700" to the catalog number for the column.

Nominal additional charge \$15.00

Example: 100mm x 4.6mm ID Ultra C18 column with Trident Integral Inlet Fitting: 9174315-700

Also order an XG-XF fitting (10mm: cat.#25026 or 20mm: 25062), see page 196.



Chromatogram Search Tool

Search by compound name, synonym, CAS # or keyword

www.restek.com/chromatograms

