



Bulk Packing Materials

Use our bulk packing materials to pack your own columns!

- Prepare your own columns in conventional or custom dimensions.
- Consistent, high-quality materials.

Restek is among the small group of column manufacturers capable of producing their own high purity silica. We offer our Ultra II®, Pinnacle® II and Pinnacle® DB silica and bonded phases in bulk. Pinnacle® II is an excellent replacement for Hypersil® silica. Pinnacle® DB is a highly base-deactivated material for analyses of basic compounds and an excellent alternative to Hypersil® BDS silica. Bonded phases from our Ultra, Allure®, and Viva lines also are available in bulk.

Our extensive QC program ensures the high quality and reproducibility of these silicas. Each lot of material is tested for mean particle size and distribution, pore diameter, surface area, and total metals analysis. You can be confident that you are getting consistent, high-quality product.

Use these materials for easy scale up to preparative scale chromatography, or for packing your own columns.



Description	min. qty.	cat.#	5-99 grams	100-499 grams	500-999 grams	≥1000 grams
3µm Ultra II Bulk Packing Materials						
Ultra II C18 Bulk Packing	5g	96043				
Ultra II C8 Bulk Packing	5g	96033				
Ultra II Silica Bulk Packing	5g	96003				
5µm Ultra II Bulk Packing Materials						
Ultra II C18 Bulk Packing	5g	96045				
Ultra II C8 Bulk Packing	5g	96035				
Ultra II Silica Bulk Packing	5g	96005				
5µm Pinnacle DB Bulk Packing Materials						
Pinnacle DB C18 Bulk Packing	5g	94145				
Pinnacle DB C8 Bulk Packing	5g	94135				
Pinnacle DB Cyano Bulk Packing	5g	94165				
Pinnacle DB Silica Bulk Packing	5g	94105				
3µm Pinnacle II Bulk Packing Materials						
Pinnacle II C8 Bulk Packing	5g	92133				
Pinnacle II C18 Bulk Packing	5g	92143				
Pinnacle II Cyano Bulk Packing	5g	92163				
Pinnacle II Phenyl Bulk Packing	5g	92153				
Pinnacle II Silica Bulk Packing	5g	92103				
5µm Pinnacle II Bulk Packing Materials						
Pinnacle II Amino Bulk Packing	5g	92175				
Pinnacle II C8 Bulk Packing	5g	92135				
Pinnacle II C18 Bulk Packing	5g	92145				
Pinnacle II Cyano Bulk Packing	5g	92165				
Pinnacle II Phenyl Bulk Packing	5g	92155				
Pinnacle II Silica Bulk Packing	5g	92105				

also available

**Restek Pack in a Box Kit:
HPLC Column Packing System**
See page 351.



Description	qty.	cat.#	1-9 bottles	10-49 bottles	50-99 bottles	≥100 bottles
5µm Ultra Bulk Packing Materials						
Ultra C1 Bulk Packing	10g/btl.	91015				
Ultra C4 Bulk Packing	10g/btl.	91025				
Ultra C8 Bulk Packing	10g/btl.	91035				
Ultra C18 Bulk Packing	10g/btl.	91745				
Ultra Amino Bulk Packing	10g/btl.	91075				
Ultra Cyano Bulk Packing	10g/btl.	91065				
Ultra Phenyl Bulk Packing	10g/btl.	91055				
Ultra Silica Bulk Packing	10g/btl.	91005				

also available

Other stationary phases and particle sizes are available; please inquire.