




Protect your column with UltraShield and UltraLine UHPLC Filters

A cost-effective way to extend the lifetime of any UHPLC column, without sacrificing UHPLC performance.



UltraShield UHPLC PreColumn Filter

- Cost-effective protection for UHPLC systems.
- Reliable way to extend column lifetime.
- Universal fit—connects easily to any brand column.
- Leak-tight to 15,000 psi (1034 bar).
- 0.5 μm titanium filter in stainless steel body with PEEK ferrule.

Specifications:

Inlet/Outlet: Female/Male 10-32
 Port Geometry: Parker ($1/16$ CPI)
 Material: Stainless Steel, PEEK ferrule
 Filter: 0.5 μm Titanium
 Pressure Rating: 15,000 psig (1054 bar)
 Wrench Flat: $5/16$ "

Description	qty.	cat.#	price
UltraShield UHPLC PreColumn Filter	ea.	24995	\$50
UltraShield UHPLC PreColumn Filter	5-pk.	24996	\$190
UltraShield UHPLC PreColumn Filter	10-pk.	24997	\$360

UltraLine UHPLC In-Line Filter

- In-line design installs easily with standard fittings.
- Cost-effective protection for UHPLC systems.
- Reliable way to extend column lifetime.
- Leak-tight to 15,000 psi (1034 bar).
- Replaceable 0.5 μm stainless steel filter in stainless steel body.

Specifications:

Inlet/Outlet: Female/Female 10-32
 Port Geometry: Parker ($1/16$ CPI)
 Material: Stainless Steel housing
 Filter: 0.5 μm Stainless Steel, 0.125" W x 0.062" T, 5 μL volume
 Pressure Rating: 15,000 psig (1054 bar)
 Wrench Flat: $3/8$ "

Description	qty.	cat.#	price
UltraLine UHPLC In-Line Filter (In-Line Assembly with Filter)	ea.	24993	\$125
UltraLine Replacement Filters	5-pk.	24994	\$60