



COLUMN INSTALLS THIS END

DI Liners for Agilent GCs** (For 0.25/0.32/0.53mm ID Columns)	Benefits/Uses	ID* OD x Length	Similar to Agilent part #	ea.	cat.#/price 5-pk.	25-pk.
Drilled Uniliner (hole near top)	trace, active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	—	21054	21055	20998
Drilled Uniliner (hole near bottom)	trace, active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	G1544-80730	20756	20771	—
Double Gooseneck Drilled Uniliner (hole near top)	trace, active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	—	20508	20509	—
Double Gooseneck Drilled Uniliner (hole near bottom)	trace, active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	G1544-80700	20954	20989	—
Drilled Cyclo-Uniliner (hole near top)	trace, dirty, high MW, active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	—	22979	22980	—
1mm Uniliner	trace, active samples, samples <1µL	1.0mm 6.3mm x 78.5mm	—	21052	21053	—
Uniliner	trace, active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	<b>NEW!</b>	22268	22269	—
Cyclo-Uniliner	trace, dirty, high MW active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	<b>NEW!</b>	22270	22271	—
Open-top Uniliner w/Wool	trace, dirty, active samples, high recovery & linearity	4.0mm 6.3mm x 78.5mm	<b>NEW!</b>	22272	22273	—

\*Nominal ID at syringe needle expulsion point.

\*\*Hole in Drilled Uniliner liner makes direct injection possible with EPC-equipped Agilent 6890 & 7890 GCs!

**i tech tip**

**Drilled Uniliner® Liners**

Use the Drilled Uniliner® liner with the hole near the bottom for semivolatiles analysis or when compounds of interest could be affected by a tailing solvent peak. Use the Drilled Uniliner® liner with the hole near the top for chlorinated pesticides analysis, aqueous injections, as well as for analysis in which the compounds of interest elute away from the solvent peak.

All liners are  
**100% deactivated**

All liners are shipped intermediate polarity (IP) deactivated unless otherwise requested.

**Sky**  
Inlet Liners

Looking for a Simple Solution?



New Sky™ liners, featuring a state-of-the-art deactivation, give you the inertness you need for accurate, reproducible trace level results.

See pages 206-211 for details.



**Siltek® Metal Inlet Liners for Agilent GCs**

- Won't crack, chip, or break like glass liners.
- Excellent response for pesticides, phenols, and other active compounds.

Liner Type (5.2mm ID x 6.3mm OD x 78.5mm)	5-pk.	25-pk.
Cyclo/Single Gooseneck	20974	20975
Single Gooseneck	21702	21703
Cyclosplitter	20726	20729
Split/Splitless w/Wool	21700	21701