




### Heater Block Retaining Nuts

for Agilent GC Split/Splitless Injectors

- Aluminum construction.
- Meet or exceed manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#	price
Heater Block Retaining Nut for Agilent 5890 GCs	19251-20620	ea.	22080	
Heater Block Retaining Nut for Agilent 6890/6850 GCs	G1544-20590	ea.	23042	



### Replacement Chemical Traps

for Agilent GCs

- Easy to install.
- Attach to same fittings as original manufacturer's equipment.
- Built-in frits retain fine particles; adsorbents remove both moisture and hydrocarbons.

Description	Similar to Agilent part #	qty.	cat.#	price
Replacement Split Vent Trap for Agilent 6890/6850 GCs	G1544-80550	ea.	22820	
Replacement Chemical Trap for Agilent 5890 GCs	05890-61260	ea.	21610	
Split Vent Line (32-inch) for Agilent GCs				
(Includes all installation hardware.)	19251-80525	2-pk.	22800	
O-Rings for Agilent Trap Fittings	5180-4181	25-pk.	22064	
Optional Split Vent Trap Assembly for Agilent 6890/6850 GCs	G1544-60610	kit	23031	
Replacement Traps (2) and O-Rings (4)	G1544-80530	kit	23032	



### Oven Temperature Sensor Assembly

for Agilent GCs

Description	Similar to Agilent part #	qty.	cat.#	price
Oven Temperature Sensor for Agilent 5890 GC	05890-61030	ea.	23040	
Oven Temperature Sensor for Agilent 6890 GC	G1530-61030	ea.	23039	



### 60 psig Backpressure Regulator Kit

Increase the versatility of your Agilent 5890 GC by replacing the existing 30 psig (207 kPa) backpressure regulator and gauge with our 60 psig (414 kPa) regulator. Enables you to use longer 60 m and 105 m columns as well as shorter 10 m columns. Includes complete instructions.

Description	Similar to Agilent part #	qty.	cat.#	price
Backpressure Regulator Kit	19246-60630	kit	20634	

also **available**

Looking for septa for Agilent GCs?

See **page 219** or visit [www.restek.com](http://www.restek.com) and search on septa.

