



22662

Replacement Diffusion Pump Fan Assemblies
for Agilent 5971/5972 MS

110/120 VAC

Description	Similar to Agilent part #	qty.	cat.#	price
Replacement Diffusion Pump Fan Assembly for Agilent 5971 MS	05971-60540	ea.	22662	
Replacement Diffusion Pump Fan Assembly for Agilent 5972 MS	05972-60540	ea.	22663	

Inland 45 is the highest quality rough pump oil you can use for your mass spectrometer.

Inland 45 Pump Oil

Recommended for most mass spectrometers.

- Ease at cold start.
- Low vapor pressure 10⁻⁷ torr.
- Nontoxic and noncorrosive.
- Compatible with buna-N, neoprene, and Viton® seals.
- Optimum vacuum pump performance.
- Lowest mass spectrometer background.
- Recommended for optimum mass spec performance.



24819

Description	Similar to Agilent part #	qty.	cat.#	price
Inland 45 Pump Oil	6040-0834	1 liter	24819	

Rough Pump Oil #19

for MSD Pumps, Oil Vacuum Pump

- Formulated from crude oil stocks known for their durability and line lubricating qualities.
- Use in Agilent 5973/5972/5971 and GCD mass spec systems, or in other manufacturers' MSD systems that require rough pump oil.
- Under average use, the oil in the foreline rough pump should be replaced every six months.



22687

Description	qty.	cat.#	price
Rough Pump Oil for MSD Pumps	1 liter	22687	

Ion Source Cleaning Powder

Use this aluminium oxide powder to clean surfaces that contact the sample or ion beam when you encounter poor sensitivity and inadequate abundances at high masses.



22685

Description	Similar to Agilent part #	qty.	cat.#	price
Ion Source Cleaning Powder	8660-0791	1 kg	22685	



27194

GC/MS Cleaning Kit

Poor sensitivity, loss of sensitivity at high masses, or high multiplier gain during an auto tune are all indicators that your mass spectrometer source may need to be cleaned. Restek has assembled all of the necessary components for cleaning and polishing your ion source.

The Restek GC/MS Cleaning Kit (cat.#s 27194, 27195) includes:

- lint free nylon gloves (small-2 pair)
- lint free nylon gloves (large-2 pair)
- lint free cotton cloth, 9 x 9 (10-pk.)
- micro mesh 4 x 6 sheet (4-pk.)
- aluminum oxide (75 gram jar)
- cotton tip applicators
- tweezers, large
- tweezers, small
- septum puller
- Dremel® tool, battery-operated (optional, 27194)
- tool kit bag

Reorder Kit (cat.# 27196) includes:

- lint free nylon gloves (small-2 pair)
- lint free nylon gloves (large-2 pair)
- lint free cotton cloth, 9 x 9 (10-pk.)
- micro mesh 4 x 6 sheet (4-pk.)

Description	qty.	cat.#	price
Mass Spec Cleaning Kit with Dremel Tool	kit	27194	
Mass Spec Cleaning Kit without Dremel Tool	kit	27195	
Mass Spec Cleaning Kit Replacement Parts Kit			
Includes: cloths, micro mesh sheets, small and large gloves	kit	27196	