



Heater Cartridge & Platinum Resistance Thermometer (PRT) Sensor
 for Agilent 5890 GCs

- Use with 5890 FID and split/splitless weldments.
- Meets or exceeds original manufacturer's performance.

Description	Similar to Agilent part #	qty.	cat.#	price
Injector/FID Heater and Injector/FID PRT Sensor	05890-61140	kit	22068	
Injector/FID Heater	19231-60620	ea.	22069	
Injector/FID PRT Sensor	19231-60660	ea.	23035	



FID Collector Housing Kits
 for Agilent GCs

- Meets or exceeds original manufacturer's performance.
- Available for Agilent 5890, 6890, or 7890 GCs.

Description	Similar to Agilent part #	qty.	cat.#	price
FID Collector Housing Kit for Agilent 5890 GCs (includes collector body, spanner nut, and silicone washer—order mount separately)	19231-20920	kit	23037	
FID Collector Housing Kit for Agilent 6890/7890 GCs (includes collector body, silicone washer, nut, and mount)	G1531-20550	kit	23045	
FID Collector Housing for Agilent 6890/7890 GCs	G1531-20740	ea.	23044	
Replacement Silicone Washers for FID Collector Housing for Agilent 5890/6890/7890 GCs	5180-4165	12-pk.	23064	



FID Collector Mount
 for Agilent GCs

- Meets or exceeds original manufacturer's performance.
- Available for Agilent 5890, 6890, or 7890 GCs.

Description	Similar to Agilent part #	qty.	cat.#	price
FID Collector Mount for Agilent 5890 GCs	19231-20930	ea.	23036	
FID Collector Mount for Agilent 6890/7890 GCs	G1531-20550	ea.	23043	



FID Base Weldment
 for Agilent GCs

- Meets or exceeds original manufacturer's performance.
- Includes brass nut.
- Available for Agilent 5890, 6850, or 6890 GCs.

Description	Similar to Agilent part #	qty.	cat.#	price
FID Base Weldment for Agilent 5890 GCs	19231-80580	ea.	23041	
FID Base Weldment, Pack Column FID for Agilent 6850/6890 GCs	G1531-80580	ea.	23052	
FID Base Weldment, Capillary Column FID for Agilent 6850/6890 GCs	G1531-80630	ea.	23053	

Note: 6890 GC connections to EPC modules are not compatible with the 7890 EPC modules.

Spanner Wrench
 for Agilent 5890/6890/6850/7890 FID Collector Assemblies

- Easily remove the nut from the FID collector without damaging the nut.
- Unique, ergonomic handle—easy to grip.
- Fits all instrument models.



Remove FID ignitor castle.



Easily loosens the nut by aligning the two pins on the bottom of the wrench with the two open slots on the nut...



...then turn counterclockwise...



...and remove.

Description	Similar to Agilent part #	qty.	cat.#	price
Spanner Wrench for Agilent 5890/6890/6850/7890 FID Collector Assembly	19231-00130	ea.	22329	