



**Replacement Cables**

for Agilent GCs, Integrators, & Autosamplers

- Considerable savings over original manufacturer's equipment.
- Instructions and wiring diagrams included.
- Manufactured with only the highest-quality components.



Connect an Agilent 5890 GC to an Agilent integrator (for second Inet integrator).

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	35900-60610	6 ft.	ea.	20650	



Connect an Agilent 5890 GC to remote-start another piece of equipment or to start the Agilent 5890 GC from that piece of equipment.

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	05890-61080	6 ft.	ea.	20657	



Connect an Agilent 5890 GC to an Agilent 3396 integrator to enable remote starts (non-Inet connection from GC to integrator).

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	03394-60560	6 ft.	ea.	20654	



Connect an Agilent 3396 integrator to a non-Agilent GC.

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	35900-60630	6 ft.	ea.	20658	




Connect an Agilent 3396 integrator to remote-start another piece of equipment or to start the Agilent 3396 integrator from that other piece of equipment.

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	03394-60540	6 ft.	ea.	20655	



Connect an Agilent 5890 GC to an Agilent 35900C Interface or to an Agilent 7673A Autosampler to enable remote starts.

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	35900-60700	6 ft.	ea.	20660	



Connect an Agilent 5890 GC to a non-Agilent integrator or standard strip chart recorder.

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	05890-60800	6 ft.	ea.	20652	



Automatic liquid sampler stop/start cable for Agilent GCs.

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	G1530-60930	7 ft.	ea.	20929	



Connect an Agilent 6890 GC to a PC 9F/9F (RS232).

Description	Similar to				
	Agilent part #	length	qty.	cat.#	price
Replacement Cable	G1530-60600	6 ft.	ea.	22694	