

Viton® O-Rings
for Shimadzu 17A, 2010, and 2014 GCs



Description	Max. temp.	Similar to Shimadzu part #	qty.	cat.#	price
Viton O-Rings for Shimadzu 17A, 2010, and 2014 GCs	250°C	036-11203-84	10-pk.	24899	

Graphite O-Rings
for Shimadzu 2010, and 2014 GCs



Description	Max. temp.	Similar to Shimadzu part #	qty.	cat.#	price
Graphite O-Rings for Split Liners	450°C	221-48393-91	5-pk.	20243	
Graphite O-Rings for Splitless Liners	450°C	221-47222-91	5-pk.	20244	



Septa (Plug)
for Shimadzu GCs

Thermolite Septa usable to 340°C inlet temp.			IceBlue Septa usable to 250°C inlet temp.			BTO Septa usable to 400°C inlet temp.		
qty.	cat.#	price	qty.	cat.#	price	qty.	cat.#	price
25-pk.	27153		50-pk.	27170		50-pk.	27118	
50-pk.	27154		100-pk.	27171		100-pk.	27119	
100-pk.	27155							



Septum Nut
for Shimadzu 17A, 2010, and 2014 GCs

- One piece design for ease of installation and removal.
- Made of clear anodized aluminum and high-quality stainless steel.

Description	Similar to Shimadzu part #	qty.	cat.#	price
Septum Nut for Shimadzu 17A, 2010, and 2014 GCs	221-41286-00, 221-44584-00	ea.	22079	



Injector Nut Kit
for Shimadzu 17A, 2010, and 2014 GCs

Includes 17A injector nut, 0.4 mm graphite ferrule, and 1/16-inch stainless steel capillary nut.

Description	qty.	cat.#	price
Injector Nut Kit for Shimadzu 17A, 2010, and 2014 GCs	kit	21895	
Siltek Treated Injector Nut Kit for Shimadzu 17A, 2010, and 2014 GCs	kit	22331	



Capillary Nut
for Shimadzu 17A, 2010, and 2014 GCs

Meets original manufacturer's performance.

Description	Similar to Shimadzu part #	qty.	cat.#	price
Capillary Nut for Shimadzu 17A, 2010, and 2014 GCs	221-41533-00	2-pk.	22688	



Restek Enhanced Capillary Nut
for Shimadzu 17A, 2010, and 2014 GCs

- Restek's design eliminates the slot, increasing lifetime and durability.
- Meets or exceeds original manufacturer's performance.

Description	Similar to Shimadzu part #	qty.	cat.#	price
Restek Enhanced Capillary Nut for Shimadzu 17A, 2010, and 2014 GCs	221-41533-00	2-pk.	20375	



5 mm Ferrules
for Shimadzu 17A GCs

- For use with packed columns.
- Graphite construction.

Description	qty.	cat.#	price
5mm Ferrules for Shimadzu 17A GCs	10-pk.	21121	



Graphite Ferrules
for Shimadzu 17A, 2010, & 2014 GCs

- Graphite 2-piece construction.
- Available in 0.4, 0.5, and 0.8 mm sizes.
- Packaged on mandrel for easy handling.

Ferrule ID	Fits Column ID	Similar to Shimadzu part #	qty.	cat.#	price
0.4mm	0.25mm and less	220-90765-00	10-pk.	24827	
0.5mm	0.32mm	221-32126-05	10-pk.	24828	
0.8mm	0.53mm	221-32126-08	10-pk.	24829	