



**Injector Wrench**

for Shimadzu 17A, 2010, and 2014 GCs

- Designed specifically for removing Shimadzu injection ports.
- High-quality stainless steel construction.

Description	Similar to Shimadzu part #	qty.	cat.#	price
Injector Wrench for Shimadzu 17A, 2010, and 2014 GCs	221-46977-00	ea.	21159	



**Open-End Wrench Set**

for use with Shimadzu 17A, 2010, and 2014 Capillary Installation Gauge

Description	qty.	cat.#	price
$\frac{1}{4}$ " x $\frac{3}{16}$ " and 10mm x 11mm Open-End Wrench Set for use with Shimadzu 17A, 2010, and 2014 Capillary Installation Gauge	ea.	22334	



**Capillary Installation Gauge**

for Shimadzu 17A, 2010, & 2014 GCs

- Seats ferrule onto column for consistent installations.
- Prevents crushed column ends.
- Made from high-quality stainless steel.
- For use with graphite ferrules.

Description	qty.	cat.#	price
Capillary Installation Gauge for Shimadzu 17A, 2010, and 2014 GCs	ea.	22333	



Everything you need in one complete kit!

22182

**Make Life Easier (MLE) Capillary Tool Kit for Shimadzu GCs**

**Includes:**

- capillary installation gauge for Shimadzu GCs
- injector wrench for Shimadzu GCs
- $\frac{1}{8}$ ",  $\frac{3}{16}$ ", and  $\frac{1}{4}$ " nylon brushes
- $\frac{1}{4}$ ",  $\frac{3}{8}$ ", and  $\frac{3}{16}$ " stainless steel wire tube brushes
- stainless steel surface brush
- 6 stainless steel jet reamers (0.25–0.65 mm OD)
- $\frac{3}{4}$ " x  $\frac{3}{16}$ " open end wrench
- $\frac{3}{8}$ " x  $\frac{7}{16}$ " open end wrench
- 6 mm x 7 mm open end wrench
- 8 mm x 10 mm open end wrench
- 16 mm x 17 mm open end wrench
- rubber-tipped slide-lock tweezers
- scoring wafers with handles
- inlet liner removal tool
- septum puller
- mini wool puller/insertor tool
- 4-inch tapered needle file
- swivel head flashlight
- mini hand drill set
- 15 cm compact steel ruler
- pocket magnifier
- high temperature string (1 meter)
- pipe cleaner (12-inch)
- cotton tip swabs (pk. of 25)

Description	qty.	cat.#	price
MLE Capillary Tool Kit for Shimadzu GCs	kit	22182	



**did you know?**  
Restek offers a full line of Parker® gas generators. Gas generators are a safe and cost-effective alternative to gas cylinders. See pages 304-308 for more information.