

Inlet Liner Seals

for Thermo Scientific TRACE PTV



Description	Max. temp.	Similar to TS part #	qty.	cat.#	price
Inlet Liner Seals for Thermo Scientific TRACE PTV	450°C	29013417	2-pk.	21392	

Graphite Sealing Ring

for Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL Instruments



Description	Max. temp.	Similar to TS part #	qty.	cat.#	price
Graphite Sealing Ring for Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL	450°C	290-334-06	ea.	21898	
	450°C	290-334-06	2-pk.	21899	

Septum Cap

for Split/Splitless Injector on Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL




Description	Similar to TS part #	qty.	cat.#	price
Septum Cap for Split/Splitless Injector on Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL	350 01 050	ea.	24971	



Septa (9.5 mm, 17 mm)

for Thermo Scientific GCs

- Preconditioned and ready to use.
- Packaged in precleaned glass jars.



Precision molded!

BTO®

Septum Diameter	25-pk./price	50-pk./price	100-pk./price
BTO Septa (usable to 400°C inlet temp.)			
9mm CenterGuide	—	27104	27105
9.5mm (3/8")	—	27106	27107
17mm CenterGuide	—	27116	27117

Due to the injection port temperatures, Restek recommends using only BTO septa in Thermo Scientific instruments.

HANDY septum size chart

Instrument	Septum Diameter (mm)
Thermo Scientific TRACE GC, 8000, 8000 TOP, GCQ w/TRACE, PTV, 8000 series	17
Finnigan (TMQ) GC 9001, GCQ, QCQ, TRACE 2000	9.5

Septa Holder Kits

for Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL

- Includes septum support and holder.
- Made from high-quality stainless steel.
- Silcosteel®-AC treated version helps with septum removal.

Description	Similar to TS part #	qty.	cat.#	price
Septa Holder for Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL	233 030 15, 350	kit	21299	
Silcosteel-AC Treated Septa Holder for Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL	233 030 15, 350	kit	24972	

Gold-Plated Liner Cap

for Split/Splitless Injector on Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL



Description	Similar to TS part #	qty.	cat.#	price
Gold-Plated Liner Cap for Split/Splitless Injector on Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL	290 042 90	ea.	22089	

Silver Seals

for Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL



Description	Similar to TS part #	qty.	cat.#	price
Silver Seal for Split/Splitless Injector	29033629	10-pk.	23057	
	29033629	20-pk.	23058	



Lab Gas Issues?

Restek has the solution!

The ProFLOW 6000 Electronic Flowmeter measures volumetric flow for gases across a range of 0.5-500 mL/min.

See **page 274** for more information.