



Adapters for Capillary Columns

on Thermo Scientific TRACE and Focus SSL

- Use same installation distance as manufacturer’s adaptors.
- Made of high-quality stainless steel.
- Siltek® treated version available for additional inertness.

Description	Similar to TS part #	qty.	cat.#	price
For use with standard 1/8" ferrules.				
Adaptor for Capillary Column on Detector Base		ea.	24916	
Adaptor for Capillary Column on Split/Splitless Injector		ea.	24917	
Siltek Treated Adaptor for Capillary Column on Split/Splitless Injector		ea.	20543	
For use with M4 ferrules.				
Adaptor for Capillary Column on Detector Base	347 25 436	2-pk.	24969	
Adaptor for Capillary Column on Split/Splitless Injector	347 05 451	ea.	24970	
Siltek Treated Adaptor for Capillary Column on Split/Splitless Injector		ea.	20544	

Nut for Terminal Fitting
for Thermo Scientific TRACE GCs



Description	Similar to TS part #	qty.	cat.#	price
Nut for Terminal Fitting for Thermo Scientific TRACE GCs	350 221 25	2-pk.	24896	

Fixing Nut for Capillary Column
for Split/Splitless Injector on Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL



Made of high-quality stainless steel.

Description	Similar to TS part #	qty.	cat.#	price
Fixing Nut for Capillary Column for Split/Splitless Injector on Thermo Scientific TRACE, 8000, 8000 TOP & Focus SSL	350 32 423	5-pk.	24973	

Graphite Ferrules
(M4 Fittings) for Thermo Scientific TRACE, 8000, 8000 TOP & Focus GCs



Ferrule ID	Fits Column ID	Similar to TS part #	Graphite 2-pk. price	Graphite 10-pk. price
0.3mm	0.10/0.15mm	—	22221	22222
0.4mm	0.18/0.28mm	29013488 (2-pk.) 29053488 (10-pk.)	20280	20281
0.5mm	0.32mm	29013487 (2-pk.) 29053487 (10-pk.)	20282	20283
0.8mm	0.45/0.53mm	29013486 (2-pk.) 29053486 (10-pk.)	20284	20285

FID Jet
for Thermo Scientific TRACE and Focus GCs



Meets or exceeds original manufacturer’s performance.

Description	Similar to TS part #	qty.	cat.#	price
FID Jet for Thermo Scientific TRACE and Focus GCs	404 043 01	ea.	23080	

Restek Electronic Leak Detector

Protect your GC and analytical column!

Don't forget to leak check after connecting all adaptors and fittings. The Restek Electronic Leak Detector will help ensure a leak-free connection!

See page 204.

www.restek.com/leakdetector