

















COLUMN INSTALLS THIS END



Splitless Liners for Varian 1075/1077 S/SL Injection Ports	Benefits/Uses:	ID* OD x Length	Similar to Varian part #	ea.	cat.#/price 5-pk.	25-pk.
 2mm Splitless	trace samples <2µL	2.0mm 6.3mm x 74mm	190010905	20721	20722	20723
 4mm Splitless	trace samples >2µL	4.0mm 6.3mm x 74mm	—	20904	20905	20906
 Double Gooseneck	trace, active samples up to 4µL	4.0mm 6.3mm x 74mm	—	20847	20848	20849
 Cyclo Double Gooseneck	trace, dirty, active samples up to 4µL	4.0mm 6.3mm x 74mm	—	20897	20898	—

Split Liners for Varian 1075/1077 S/SL Injection Ports	Benefits/Uses:	ID* OD x Length	Similar to Varian part #	ea.	cat.#/price 5-pk.	25-pk.
 1mm Split	purge & trap inlet splitting or samples <1µL	1.0mm 6.3mm x 72mm	—	20970	20971	—
 Splitter w/Wool	universal, use with rapid autosamplers	4.0mm 6.3mm x 72mm	190010901	20792	20793	20794
 Laminar Cup Splitter	high MW compounds	4.0mm 6.3mm x 72mm	190010902	20803	20804	—
 Cup Splitter	high & low MW compounds	4.0mm 6.3mm x 72mm	—	20724	20725	—
 Cyclo-splitter	dirty samples, many injections before cleaning required	4.0mm 6.3mm x 72mm	—	20727 \$55	20728	—
 Frit Splitter	dirty samples, non-active compounds	4.0mm 6.3mm x 72mm	190010903	20715	20716	20717
 Baffle Splitter	compounds with close boiling points	4.0mm 6.3mm x 72mm	190010904	20718	20719	20720
 Split Precision Liner w/Wool	dirty samples, active samples	4.0mm 6.3mm x 72mm	—	21030	21031	—

DI Liners for Varian 1075/1077 S/SL Injection Ports (0.32/0.53mm ID)	Benefits/Uses:	ID* OD x Length	Similar to Varian part #	ea.	cat.#/price 5-pk.	25-pk.
 Uniliner	trace, active samples, high recovery & linearity	4.0mm 6.3mm x 72mm	—	20345	20346	—
 Cyclo-Uniliner	trace, dirty, high MW, active samples, linearity	4.0mm 6.3mm x 72mm	—	20347	20348	—
 Open-top Uniliner w/Wool	trace, dirty, active samples, high recovery & linearity	4.0mm 6.3mm x 72mm	—	20845	20846	—

SPME Liners for Varian 1075/1077 S/SL Injection Ports	Benefits/Uses:	ID* OD x Length	Similar to Varian part #	ea.	cat.#/price 5-pk.	25-pk.
 SPME Liner	ideal for low-volume SPME applications	0.75mm 6.3mm x 74mm	—	21112	21113	—

*Nominal ID at syringe needle expulsion point.

All liners are
100%
deactivated

All liners are shipped intermediate polarity (IP) deactivated unless otherwise requested.