

Injection Port Conversion Chart

Instrument	Uniliner Liner Set-Up		Vu-Tight Fitting	
	Uniliner Liner: cat.# 20301, 20305, 20308, 20309, 20315, 20316, 20319, 20320, 20902, & 20903	Liner Adaptor: cat.# 20310	Vu-Tight Fitting Kit: cat.# 20504 Liner: cat.# 20342, 20343, 20344, 20787, & 20788	Vu-Tight Fitting Kit: cat.# 20504
Agilent GCs (1/4-inch injection ports) Models: 5700, 5710, 5711, 5712, 5830, 5840, 5880, 5890, 6850, 6890	✓	✓	✓	✓
Varian GCs (1/4-inch injection ports) Models: 1200, 1400, 2100, 2400, 3300-3700, 4400, 4600, 6000	✓	✓	✓	✓
Tracor GCs (1/4-inch injection ports) Models: 540, 550, 560, 565, 570	✓	✓	✓	✓
Packard Becker GCs (1/4-inch injection ports) Models: 427, 428, 429, 430, 433, 436, 437, 438	✓	✓	✓	✓
Gow-Mac GCs (1/4-inch injection ports) Models: 69-750, 69-550	✓	✓	✓	✓
HNU GCs (1/4-inch injection ports) Models: 301, 401, 421	✓	✓	✓	✓
PerkinElmer GCs (1/4-inch injection ports) Models: Sigma, 1B-4B, 300, 2001, 2100	✓	✓	✓	✓
PerkinElmer Auto SYS*	✓			

\*Does not require Uniliner liner adaptor.

**Sky**  
Inlet Liners

Introducing **Sky™** Liners

RESTEK

See page 206 for more details.