

Restek Electronic Leak Detector

Don't let a small leak turn into a costly repair—protect your instrument and analytical column by using a Restek Leak Detector.

Features & Benefits include:

- Optimized sample flow path.
- New ergonomic, hand-held design.
- Rugged side grips for added durability.
- Handy probe storage for cleanliness and convenience.
- Longer lasting battery, up to 6 hours of continuous use.
- Automatic shut-off.
- A convenient carrying and storage case.
- Easy to clean probe assembly.
- A universal charger set (US, European, UK, and Australian plugs included).

Backed by a 1-year warranty, the new Restek Leak Detector sets an industry standard for performance and affordability in hand-held leak detectors.

Leak Detector Facts

Detectable gases:	helium, nitrogen, argon, carbon dioxide, hydrogen
Battery:	Rechargeable Ni-MH internal battery pack (6 hours normal operation)
Operating Temperature Range:	32°-120°F (0°-48°C)
Humidity Range:	0-97%
Warranty:	one year
Certifications:	CE, Ex, Japan
Compliance:	WEEE, RoHS

Limits of Detection

These gases can be detected with the Restek Electronic Leak Detector at the following leak rates.

Gas	Minimum Detectable Leak Rate (atm cc/sec.)	Indicating LED Light Color
Helium	1.0 X 10 ⁻⁵	Red
Hydrogen*	1.0 X 10 ⁻⁵	Red
Nitrogen	1.4 X 10 ⁻³	Yellow
Argon	1.0 X 10 ⁻⁴	Yellow
Carbon Dioxide	1.0 X 10 ⁻⁴	Yellow



22839



22839

Hard-sided carrying case included with purchase of unit.



22657

Optional soft-side storage case is ideal for storing your leak detector in smaller spaces, such as your tool box.

Description	qty.	cat.#	price
Leak Detector with Hard-Sided Carrying Case and Universal Charger Set (US, UK, European, Australian)	ea.	22839	
Leak Detector Routine Maintenance Review**	ea.	22839-R	
Soft-Side Storage Case	ea.	22657	
Small Probe Adaptor	ea.	22658	

*Caution: The Restek Electronic Leak Detector is designed to detect trace amounts of hydrogen in a noncombustible environment. It is NOT designed for determining leaks in a combustible environment. A combustible gas detector should be used for determining combustible gas leaks under any condition. The Restek Electronic Leak Detector may be used for determining trace amounts of hydrogen in a GC environment only.

**Routine maintenance includes inspection of the probe tip, internal/external tubing and a battery replacement.



Checking for leaks maintains the integrity of a GC system.

i tech tip

To keep your Leak Detector performing at optimal levels, send it in every two years for a Routine Maintenance Review. Our engineers will inspect the probe tip, internal/external tubing and replace the battery. This preventive maintenance plan will help ensure your leak detector operates properly for years to come.

for more info

How much sensitivity is needed in a leak detector?
www.restek.com/cat001



Verify hard-to-reach leaks using the small probe adaptor.

22658

i tech tip

Avoid using liquid leak detectors on a GC! Liquids can be drawn into the system.

