



22185



22186



22182



22184



22183



MLE Tool Kits conveniently provide the tools that make it easier to install and maintain capillary columns at a discounted price compared to buying the tools individually!

### Make Life Easier (MLE) Capillary Tool Kits

Everything you need in one complete kit!

**All kits include these components:**

- 1/8", 3/16", 1/4" nylon brushes
- 1/4", 3/8", 3/16" stainless steel wire tube brushes
- stainless steel surface brush
- 6 stainless steel jet reamers (0.25–0.65 mm OD)
- 1/4" x 5/16" open end wrench
- 3/8" x 7/16" open end wrench
- rubber-tipped slide-lock tweezers
- scoring wafers with handles
- inlet liner removal tool
- septum puller
- mini wool puller/insertor tool
- 4-inch tapered needle file
- swivel head flashlight
- mini hand drill set
- 15 cm compact steel ruler
- pocket magnifier
- high temperature string (1 meter)
- pipe cleaner (12-inch)
- cotton tip swabs (pk. of 25)

**MLE Capillary Tool Kit for Agilent GCs (cat.# 22186)**

**also includes:**

- capillary installation gauge for Agilent GCs
- injector wrench for Agilent GCs
- septum nut removal tool
- 7/16" x 1/2" open end wrench
- 1/2" x 9/16" open end wrench

**MLE Capillary Tool Kit for PerkinElmer GCs**

**(cat.# 22185) also includes:**

- 7/16" x 1/2" open end wrench
- 1/2" x 9/16" open end wrench

**MLE Capillary Tool Kit for Shimadzu GCs**

**(cat.# 22182) also includes:**

- capillary installation gauge for Shimadzu GCs
- injector wrench for Shimadzu GCs
- 6 mm x 7 mm open end wrench
- 8 mm x 10 mm open end wrench
- 16 mm x 17 mm open end wrench

**MLE Capillary Tool Kit for Thermo Scientific GCs**

**(cat.# 22183) also includes:**

- capillary installation gauge for Thermo Scientific GCs
- liner cap removing tool for Thermo Scientific GCs
- 6 mm x 7 mm open end wrench
- 8 mm x 10 mm open end wrench
- 16 mm x 17 mm open end wrench

**MLE Capillary Tool Kit for Varian GCs (cat.# 22184)**

**also includes:**

- capillary installation gauge for Varian GCs
- 7/16" x 1/2" open end wrench
- 1/2" x 9/16" open end wrench

Description	qty.	cat.#	price
MLE Capillary Tool Kit for Agilent GCs	kit	22186	
MLE Capillary Tool Kit for PerkinElmer GCs	kit	22185	
MLE Capillary Tool Kit for Shimadzu GCs	kit	22182	
MLE Capillary Tool Kit for Thermo Scientific GCs	kit	22183	
MLE Capillary Tool Kit for Varian GCs	kit	22184	