



Restek Filter Base Plates

- Standard base plate fittings are 1/8". To adapt to 1/4", order 1/8" to 1/4" tube-end unions.
- End fittings available in brass or stainless steel.
- Base plates fit all Super-Clean gas filters listed.

| Description | qty. | Brass | | Stainless Steel | |
|------------------------------------|------|-------|-------|-----------------|-------------|
| | | cat.# | price | qty. | cat.# price |
| Filter Base Plate, Single-Position | ea. | 22025 | \$231 | ea. | 22344 |
| Filter Base Plate, 2-Position | ea. | 22026 | \$425 | ea. | 22345 |
| Filter Base Plate, 3-Position | ea. | 22027 | \$609 | ea. | 22346 |

Wall Mounting Bracket

Base plates can be mounted by using screws and the mounting holes on the base plate, or by using this optional wall mounting bracket.



| Description | qty. | cat.# | price |
|---|------|-------|-------|
| Wall Mounting Bracket for Super-Clean Base Plates | ea. | 21984 | |

Replacement O-Rings for Cartridge Base Plates

Pack includes 10 large O-rings and 10 small O-rings.



| Description | qty. | cat.# | price |
|---|--------|-------|-------|
| Replacement O-Rings for Cartridge Base Plates | 20-pk. | 22023 | |

1/8-Inch to 1/4-Inch Tube-End Unions

To adapt 1/8" Super-Clean base plate fittings to 1/4", use 1/8" to 1/4" tube-end unions.



| Description | qty. | Brass | | Stainless Steel | |
|-------------------------------------|-------|-------|-------|-----------------|-------------|
| | | cat.# | price | qty. | cat.# price |
| Tube-End Reducer, 1/8" tube to 1/4" | 5-pk. | 21833 | \$31 | 2-pk. | 21933 |

Restek Super-Clean Gas Trapping System for LC/MS

A Super-Clean quick-change cartridge system efficiently removes hydrocarbons from nitrogen!

- Changing filters is quick and easy.
- Up to 20 L of hydrocarbon-free nitrogen per minute.
- Filters connected in parallel to handle high flows of LC/MS.



20 L of purified nitrogen per minute!

Super-Clean Gas Filters provide high-purity outlet gas

| | |
|-----------------------|------------------------|
| Type of filter: | Hydrocarbon (charcoal) |
| Max. Flow: | 20 L/min. |
| Outlet Gas Quality %: | 99.9999% |
| Maximum Pressure: | 11 bar/159 psi |
| Estimated Lifetime: | 3 to 6 months |

| Description | qty. | cat.# | price |
|---|-------|-------|-------|
| Super-Clean Gas-Trapping System (2-position base plate, 2 charcoal filters) | ea. | 22062 | |
| Replacement Hydrocarbon (Charcoal) Filters | 2-pk. | 22061 | |

also **available**

Looking for a nitrogen generator for your LC/MS?

Restek offers a full line of Parker LC/MS generators.

See **page 307**.



Catch the Buzz

Sign up for Restek's e-newsletter, *The Buzz*

www.restek.com/buzz