

Tubing Tools



1/16-Inch Tubing Cutter

20192

- Produces square, smooth cuts in 1/16-inch tubing.
- Eliminates tubing distortion.
- Replaceable cutting wheel.

Description	qty.	cat.#	price
1/16" Tubing Cutter	ea.	20192	
Replacement Cutting Wheels	3-pk.	20185	



Dimensions: 8" x 6 1/4" x 4 1/4"
(20.3 x 15.9 x 10.8 cm)
Weight: 11 lbs. (5.0 kg)

SSI TC-20 Tube Cutting Machine

- Cuts 1/16", 1/8", or 1/4" tubing with inside diameter as small as 0.008".
- Electrically operated bench-top model.
- Handy dressing tool on the swing arm removes burrs and reams tubing.
- Voltage selectable 110–120/220–240 volts, 50–60 Hz.*

Description	qty.	cat.#	price
SSI Tubing Cutter Machine	ea.	23029	
SSI Replacement Cutting Wheels	3-pk.	23030	

*Unit shipped set for 110–120 operating voltage. Switch to 220–240 volts by using alternate fuse and power cord (included).



Ridgid® Tubing Cutter

- Excellent for cutting 1/8- or 1/4-inch metal tubing.*
- Compact size is ideal for tight spaces.
- Replaceable cutting wheel.

Description	qty.	cat.#	price
Ridgid Tubing Cutter for 1/8" or 1/4" metal tubing	ea.	23011	
Replacement Cutting Wheels	2-pk.	23012	

*Not for use with stainless steel tubing.



Ridgid® Heavy-Duty Tubing Cutter

- Specifically designed to cut 3/16-1 1/8" stainless steel tubing.
- Large knob and 6 individual bearings for more control with less turning.
- Convenient fold-away reamer.

Description	qty.	cat.#	price
Ridgid Heavy-Duty Tubing Cutter	ea.	22621	



1/16-Inch Tubing Cutting Pliers

- Ideal for cutting 1/16-inch tubing.
- Cuts quickly, reducing distortion.
- Cuts clean, eliminating need for deburring.


Description	qty.	cat.#	price
1/16" Tubing Cutting Pliers	ea.	20193	

Restek Tubing Scorer

for MXT® Columns

- Makes perfect cuts every time.
- Easy to use.
- Leaves column entrance perfectly round.

Metal MXT® columns are easy to cut. Scoring wafers can be used, but may leave the column end irregularly shaped. The Restek tubing scorer is designed to make a perfect cut every time, leaving the column entrance perfectly round.



Description	qty.	cat.#	price
Restek Tubing Scorer for MXT Columns (0.25-0.53mm ID & 0.5-0.8mm OD)	ea.	20523	
Replacement Scoring Wheel	ea.	20522	