



25322

Kits are convenient and economical, compared to ordering individual items!



25097

Survival Kit for HPLC, PEEK™

The PEEK™ Survival Kit is an invaluable parts kit that contains tubing, fittings, and tools essential for setting up and maintaining your HPLC system: PEEK™ tubing, connectors, and elbows, Teflon® tubing, a tubing cutter and extra blades, a ValvTool wrench, open-end wrenches, and more.

For start-up and maintenance in all HPLC systems.

Description	qty.	cat.#	price
Survival Kit for HPLC, PEEK	kit	25322	

Contents:

Description	qty.	cat.#	price
PEEK Column Connector (beige, round body)	10-pk.	25015	
PEEK Tubing, 1/16" OD x 0.005" ID x 3m (red stripe)	ea.	25065	
PEEK Tubing, 1/16" OD x 0.007" ID x 3m (yellow stripe)	ea.	25066	
PEEK Tubing, 1/16" OD x 0.010" ID x 3m (blue stripe)	ea.	25067	
PEEK Tubing Elbow, 90°	5-pk.	25308	
PEEK Tubing Elbow, 180°	5-pk.	25309	
Teflon Tubing, 1/8" OD x 0.063" (1.6mm) ID x 3m	3m	25306	
Teflon Tubing, 1/8" OD x 0.094" (2.4mm) ID x 3m	3m	25307	
Tubing Clip	5-pk.	25310	
ValvTool Wrench	ea.	25321	
Open-End Wrenches, 1/4" x 5/16"	2-pk.*	20110	
Clean-Cut Tubing Cutter	ea.	25069	
Replacement Blade for Clean-Cut Cutter	ea.	25070	
PEEK Union Connector 1/16"	2-pk.	25323	
Mobile Phase Sparge Filter: 2µm	ea.	25311	
Inlet Filter: 10µm	ea.	25312	

*Kit contains 1 wrench, replacement (cat.# 20110) is a 2-pk.

HPLC PEEK™ Performance Kit

This HPLC PEEK™ Performance Kit is a useful complement to your HPLC system. Includes:

- Hub-Cap Adapter w/Cap (1)
- Tubing Clips (3)
- Clean Cut Tubing Cutter (1)
- PEEK Finger Tight Fittings (5)
- PEEK Column Plugs (5)
- PEEK 90° Tubing Elbows (2)
- PEEK 180° Tubing Elbows (2)
- PEEK Column Connectors (2)



Description	qty.	cat.#	price
HPLC PEEK Performance Kit	kit	26394	

Survival Kit for HPLC, Stainless Steel

The stainless steel survival kit contains a wide range of tubing, fittings, and tools necessary to set up and maintain your HPLC system: a selection of lengths and IDs of 1/16-inch tubing, plus nuts, ferrules, a ValvTool wrench, and a zero-dead-volume union.

Description	qty.	cat.#	price
Survival Kit for HPLC, Stainless Steel	kit	25097	

Contents:

Description	color	qty.	cat.#	price
HPLC Capillary Tubing, SS 1/16" x 0.005" x 5cm	red	3-pk.	25240	
HPLC Capillary Tubing, SS, 1/16" x 0.005" x 10cm	red	3-pk.	25241	
HPLC Capillary Tubing, SS, 1/16" x 0.005" x 20cm	red	3-pk.	25242	
HPLC Capillary Tubing, SS, 1/16" x 0.005" x 30cm	red	3-pk.	25243	
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 5cm	gray	3-pk.	25244	
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 10cm	gray	3-pk.	25245	
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 20cm	gray	3-pk.	25246	
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 30cm	gray	3-pk.	25247	
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 5cm	blue	3-pk.	25248	
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 10cm	blue	3-pk.	25249	
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 20cm	blue	3-pk.	25250	
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 30cm	blue	3-pk.	25251	
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 5cm	yellow	3-pk.	25252	
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 10cm	yellow	3-pk.	25253	
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 20cm	yellow	3-pk.	25254	
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 30cm	yellow	3-pk.	25255	
1/16" Rheodyne Style Nut		10-pk.	25095	
1/16" Rheodyne Style Ferrule		10-pk.	25096	
ValvTool Wrench		ea.	25321	
Ferrules, 1/16" Stainless Steel		10-pk.	26548	
Nuts, 1/16" Stainless Steel		10-pk.	20287	
Zero-Dead-Volume Internal Union		ea.	20148	