



Now available for GL-45 bottles!



Bottles not included

**Waste Overflow Indicator**

- Avoid messy pooling around mobile phase waste containers.
- Audible alarm instantly alerts user, preventing overflow.
- Compact, battery operated unit.
- Available for 4-liter and GL-45 solvent bottles.

The Restek Waste Overflow Indicator will help keep your mobile phase waste where it belongs—in the waste container! Compact, battery operated unit fits securely on solvent bottles and accommodates two waste streams. An audible alarm is given as the solvent waste container approaches capacity, giving you time to empty or change the container. Another innovative design from Restek!

Description	qty.	cat.#	price
Waste Overflow Indicator for HPLC Systems, 4 Liter	ea.	26543	
Waste Overflow Indicator for HPLC Systems, GL-45	ea.	26550	
Replacement AA Battery for the Waste Overflow Indicator	ea.	26544	
Replacement AA Batteries for the Waste Overflow Indicator	3-pk.	26545	

**Opti-Cap® Bottle Top**

The most economical way to helium-sparge and deliver HPLC mobile phases.

The Opti-Cap® top fits all standard GL-45 bottles and has two 1/8-inch holes and one 1/16-inch hole for tubing. All three openings are designed for threaded PEEK™ plugs



25303

Bottle not included.

Description	qty.	cat.#	price
Opti-Cap (Cap and PEEK Plug)	ea.	25300	
Opti-Cap Kit (Opti-Cap, 3 meters of tubing, sparging filters)	kit	25301	
Opti-Cap Kit with 1L Bottle	kit	25302	
Opti-Cap Kit with 2L Bottle	kit	25303	
Related items and replacement parts	qty.	cat.#	price
Mobile Phase Sparge Filter: 2µm, stainless steel	ea.	25311	
Mobile Phase Inlet Filter: 10µm	ea.	25312	
Teflon Tubing, 1/8" OD x 0.094" (2.4mm) ID x 3m	3m	25307	
Teflon Tubing, 1/8" OD x 0.063" (1.6mm) ID x 3m	3m	25306	
PEEK Plug, 1/4"-28 threads	3-pk.	25319	
1L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25304	
2L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25305	

**Eco-Cap Bottle Top, PTFE**

Fits all standard GL-45 bottles; has two 1/8-inch holes and one 1/16-inch hole for tubing. Holes are not threaded.



25395

Description	qty.	cat.#	price
Eco-Cap Bottle Top, includes 1 Male Luer Plug	ea.	25395	



27197

26551

**Opti-Cap® Adapters**

Allow the use of the Hub-Cap with GL-45 solvent bottles.

Description	qty.	cat.#	price
Opti-Cap Adapter	ea.	27197	
Opti-Cap Adapter Multi-pack	3-pk.	27198	
Opti-Cap Adapter w/Hub-Cap	kit	26551	