

Preassembled Vial Kits and Limited Volume Inserts for 2.0 mL Big Mouth Step Design Screw-Thread Vials, 12 x 32 mm, 10 mm/425 Thread

2.0 mL Big Mouth Step Vial Kits w/Black Caps (preassembled)
Caps attached to vials.

Description	100-pk.	1000-pk.
Clear 2.0mL Vial, Black Cap, PTFE/Silicone/Red PTFE Septa, 0.040"	24675	24676
Clear 2.0mL Vial, Black Cap, Red PTFE/Silicone Septa, 0.065"	24642	24643
Amber 2.0mL Vial, Black Cap, Red PTFE/Silicone Septa, 0.065"	24644	24645

Inserts
for 2.0 mL Big Mouth Step Design Screw-Thread Vials

Description	Volume	Material	100-pk.	1000-pk.
Big Mouth Insert w/Bottom Spring	250µL	Glass	21776	21777
Insert, Flat Bottom	350µL	Glass	21780	24517
Insert, Flat Bottom w/ID Ring	350µL	Glass	24692	24693
Big Mouth Insert w/Glass Flange (Step Design)*	250µL	Glass	24516	21779

Polypropylene inserts available on request (1,000-packs only).
*Not to be used with 9mm screw-thread vials.



WISP Vials

4.0 mL WISP 48 Crimp-Top 13 mm Vials, 15 x 45 mm, Snap Seal

4.0 mL WISP 48 Crimp-Top Vial (vial only)

Description	Volume	Color	100-pk.	1000-pk.
WISP 48 Snap Seal Vial w/White Graduated Marking Spot	4.0mL	Clear	24658	24659



24658

Aluminum Seals w/Septa
for 4.0 mL WISP Crimp-Top Vials

Seal Color	Septa Material	100-pk.	1000-pk.
Silver	PTFE/Natural Rubber	21753	21754
Silver	PTFE/Silicone	21755	21756



also **available**
See **page 362** for inserts.

instrument reference

For 4.0 mL WISP 48 Style Vials

The following instruments are compatible with 4.0 mL WISP 48 style vials. For further information or questions, please contact your Restek sales representative or technical service.

Manufacturer	Instrument/Model #	Manufacturer	Instrument/Model #
Bruker	LC 51	TosoHaas	TSK-6080, AS-8010, AS-8020
Hitachi	L-7200, L7250	Unicam	LC-XP
Phillips	LC-XP	Waters	48-pos. M700, 710, 717+
Shimadzu	AOC-14/1400, AOC-17, SIL-2AS		