



22701

2.0 mL Versa Vials™, Glass, 12 x 32 mm

- Increased neck length for accommodating fluctuations in robotic arm settings.
- Eliminate cracked necks, breakage, and dangerous sharp edges, using soft Versa Vial PTFE/silicone plugs.
- Largest target area available, ideal for sampling and pipet use.
- Vials available with or without graduated marking spots.
- Patent pending.

Description	Volume	Color	Material	100-pk.	1000-pk.
Versa Vial	2.0mL	Clear	Glass	22701	22702
Versa Vial	2.0mL	Amber	Glass	22697	22704
Versa Vial w/Graduated Marking Spot	2.0mL	Clear	Glass	22703	22704



22705

1.5 mL Versa Vials™, Polypropylene, 12 x 32 mm

Description	Volume	Color	Material	100-pk.	1000-pk.
Versa Vial	1.5mL	Clear	Polypropylene	22705	22706



22707

Versa Vial™ Inserts

Description	Volume	Material	100-pk.	1000-pk.
Versa Vial Insert w/Glass Flange	250µL	Glass	22707	—
Versa Vial Insert w/Flange	250µL	Polypropylene	22709	22710
Versa Vial Insert, Flat Bottom	250µL	Glass	22711	22712



22715

Versa Vial™ Plugs (12 mm)

Description	Color	Plug Material	100-pk.	1000-pk.
Versa Vial Plug	Green	Polyethylene	22713	22714
Versa Vial Plug	Grey	Chlorobutyl/Siliconized	22715	22716
Versa Vial Plug	White	PTFE/Silicone	22717	22718
Versa Vial Plug	White	PTFE/Silicone Plug w/Slit*	22719	22720

*Eliminates needle damage to plug.



21798

26122

Pre-Cleaned Volatile Organic Analyte (VOA) Sampling Vials

- Container, liner, and closure cleaned and assembled.
- Clear or amber.
- Open top caps.
- Teflon® faced 0.125" silicone septa.
- Each case lot numbered.

Description	Volume	Color	Material	Screw-Thread Size	qty.	cat.#	price
Precleaned VOA Vials	20mL	Clear	Glass	24mm/400	72-pk.	21798	
Precleaned VOA Vials	20mL	Amber	Glass	24mm/400	72-pk.	21799	
Precleaned VOA Vials	40mL	Clear	Glass	24mm/400	72-pk.	21796	
Precleaned VOA Vials	40mL	Amber	Glass	24mm/400	72-pk.	21797	
Collection Vials for ASE 200	60mL	Clear	Glass	24mm/400	72-pk.	26121	
Collection Vials for ASE 200	60mL	Amber	Glass	24mm/400	72-pk.	26122	
Replacement Septa for 20, 40 & 60mL Vials and 250mL Bottles, 24mm x 0.125"			Teflon-lined silicone		100-pk.	24694	