

**Gas-Tight Teflon®-Tipped Syringes for Agilent 7673, 7683, and 6850 Autosamplers**

**HAMILTON**

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		Restek		price
					Model	cat.#	qty.	cat.#	
10µL	ASN	23s	1.71"	Agilent	1701	80080	ea.	24894	
10µL	ASRN	23s	1.71"	Agilent	1701	80087	ea.	21224	
10µL	ASRN	26s	1.71"	Agilent	1701	80088	ea.	21228	
10µL	ASRN	23s-26s	1.71"	Agilent	1701	80089	ea.	21229	

**SGE**

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		Restek		price
					Model	cat.#	qty.	cat.#	
10µL	F	23	42mm	Cone	10F-HP-GT-0.63	002812	ea.	24789	
10µL	R	23	42mm	Cone	10R-HP-GT-0.63	002818	ea.	24797	
10µL	R	23-26	42mm	Cone	10R-HP-GT-0.63/0.47	002829	ea.	21220	

**Large Volume Syringes for Agilent 7673, 7683, and 6850 Autosamplers**

**SGE**

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		Restek		price
					Model	cat.#	qty.	cat.#	
25µL	R	23	42mm	Cone	25R-HP-0.63	003665	ea.	24798	
50µL	R	23	42mm	Cone	50R-HP-0.63	004665	ea.	24799	
100µL	R	23	42mm	Cone	100R-HP-0.63	005665	ea.	24800	
250µL	R	23	42mm	Cone	250R-HP-0.63	006665	ea.	24801	



24894

**Replacement Needles for Agilent 7673, 7683, and 6850 Autosampler Syringes**

**HAMILTON**

Syringe Volume	Needle Gauge	Needle Length	Point Style	Hamilton		Restek		price
				cat.#	qty.	cat.#	qty.	
5-10µL	23s	1.71"	Agilent	7786-01	6-pk.	20164		
5-10µL	26s	1.71"	Agilent	7786-02	6-pk.	20165		
5-10µL	23s-26s	1.71"	Agilent	7785-01	6-pk.	20166		

**SGE**

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE		Restek		price
				cat.#	qty.	cat.#	qty.	
0.5µL	23	42mm	Cone	033715	ea.	24821		
0.5µL	26	42mm	Cone	033708	ea.	24820		
5µL	23	42mm	Cone	036720	2-pk.	24823		
10µL	23	42mm	Cone	037717	2-pk.	24825		
10µL	26	42mm	Cone	037715	2-pk.	24824		
25-250µL	23	42mm	Cone	038717	2-pk.	24826		

Do not reuse TFE washers.

also **available**

**Looking for an Agilent Syringe Cross-Reference?**

See page 372.

