



**Replacement Plunger Assembly, Teflon®-Tipped for Agilent 7673, 7683, and 6850 Autosampler Syringes**  
Teflon®-Tipped


**SGE**

Syringe Volume	Needle Term.	SGE cat.#	qty.	Restek cat.#	price
10µL	F	031808	2-pk.	21218	
10µL	R	031809	2-pk.	21284	

**Needle Gauge for Agilent 7673, 7683, 6850 Syringes**


**Single Gauge Needle (23s)**

- The most popular gauge for Agilent 7673.
- Best for Merlin Microseal septum and standard septum-equipped GCs.
- Packed column injection ports.
- Split/splitless injection ports.




**Single Gauge Needle (26s)**

- On-column injection ports.
- Split/splitless injection ports.



**23s-26s—Dual Gauge (tapered)**

- Durability of a 23s gauge needle.
- Ability of a 26s gauge needle to perform split/splitless and on-column injections.



**Cross-reference**

**For Agilent Syringes**

Listed below are some Hamilton and SGE syringes, cross-referenced to Agilent part numbers for your convenience. For further information or questions, contact your Restek sales representative or technical service.

Similar to Agilent Part #	Restek Part #	Description	Page
5181-1267	21212, 24596	10µL Standard Syringe, Fixed Needle, 23-26 Gauge	370
5181-1273	21210, 24593	5µL Standard Syringe, Fixed Needle, 23-26 Gauge	370
5181-3319	20166	Needle for 10µL Standard Syringe, 23-26 Gauge	371
5181-3321	21213, 24598	10µL Standard Syringe, Removable Needle, 23-26 Gauge	370
5181-3356	21218	Replacement Teflon-Tipped Plunger for 10µL Fixed Needle Syringe	372
5181-3358	24919	Teflon-Tipped Plunger for 10µL Standard Syringe, Removable Needle	382
5181-3360	21215, 24600	10µL Standard Syringe, Fixed Needle, 23-26 Gauge (6-pk.)	370
5181-3365	21218	Teflon-Tipped Plunger for 10µL Standard Syringe, Fixed Needle	372
5181-8809	24894, 24789	10µL Gas-tight Syringe, Fixed Needle, 23 Gauge	371
5181-8811	20164, 24825	Needle for 10µL Gas-tight Syringe, 23 Gauge	371
5181-8813	21224, 24797	10µL Gas-tight Syringe, Removable Needle, 23 Gauge	371
5182-0830	24821	Needle for 0.5µL Standard Syringe, 23 Gauge	371
5182-0832	20166	Needle for 5-10µL Autosampler Syringe, 23-26 Gauge	371
5182-0834	20172, 24793	5µL Standard Syringe, Removable Needle, 23 Gauge	370
5182-0835	21223, 21211	5µL Standard Syringe, Removable Needle, 23-26 Gauge	370
5182-0875	20170, 24783	5µL Standard Syringe, Fixed Needle, 23 Gauge (6-pk.)	370
9301-0713	20167, 24785	10µL Standard Syringe, Fixed Needle, 23 Gauge	370
9301-0714	24595, 24784	10µL Standard Syringe, Fixed Needle, 26 Gauge	370
9301-0725	20169, 24787	10µL Standard Syringe, Fixed Needle, 23 Gauge (6-pk.)	370
9301-0891	24592, 24780	5µL Standard Syringe, Fixed Needle, 26 Gauge	370
9301-0892	20168, 24781	5µL Standard Syringe, Fixed Needle, 23 Gauge	370