



Autosampler Syringes for Varian 8000 Series GCs

- Hamilton and SGE syringes are designed and tested to meet critical autosampler performance needs.
- Needle point styles are designed to withstand multiple, rapid injections through a septum.

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	Hamilton cat.#	Restek cat.#	Restek price
Varian 8100/8200 Syringes								
10µL	RN	26s	1.99"	5	701 Varian	202880	24956	ea.

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	SGE cat.#	Restek cat.#	Restek price
Varian 8000 Syringes								
10µL	R	25	50mm	Dome	10R-VA80-II	002926	24850	ea.
Varian 8035, 8100, & 8200 Syringes								
10µL*	F	25	53mm	S/Hole	10F-VA8X-11	002923	24851	ea.
10µL*	R	25	53mm	S/Hole	10R-VA8X-II	002924	24852	ea.
100µL	R	25	53mm	S/Hole	100R-VA8X	005921	24853	ea.

Varian CP 8400, CP8410 Syringes

10µL	F	26	50mm	Cone	10F-C/F-5/0.47C	002980	24922	ea.
10µL	F	26	50mm	Bevel	10F-VA8400-5/0.47	002950	21202	ea.

*Gas-tight syringe.

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) Teflon® Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip



Replacement Needles for Varian 8000 Series Autosampler Syringes

SGE

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	Restek cat.#	Restek price
10µL	23	53mm	S/Hole	037779	24857	2-pk.
10µL*	25	53mm	S/Hole	037780	24858	ea.
10µL**	25	53mm	S/Hole	037777	24854	ea.
10µL	25	53mm	2	037776	24855	2-pk.

*straight needle shaft, 0.20mm ID needle, for use with viscous samples

**tapered needle shaft, 0.12mm ID needle



also **available**

Looking for Restek Genuine Replacement Parts for Varian Instruments?
See pages 259-264.