

Teflon® Tip, Gas-Tight Syringes

- Suitable for gases or liquids, for maximum inertness.
- High accuracy of dispensed volumes.
- Interchangeable barrels, plungers, and tips extend performance and increase cost-effectiveness.

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		Restek	price
					Model	cat.#		
10µL	N	26s	2"/51mm	2	1701	80000	ea.	24557
10µL	RN	26s	2"/51mm	2	1701	80030	ea.	24558
25µL	N	22s	2"/51mm	2	1702	80200	ea.	24559
25µL	RN	22s	2"/51mm	2	1702	80230	ea.	24560
50µL	N	22s	2"/51mm	2	1705	80900	ea.	24561
50µL	RN	22s	2"/51mm	2	1705	80930	ea.	24562
100µL	N	22s	2"/51mm	2	1710	81000	ea.	24563
100µL	RN	22s	2"/51mm	2	1710	81030	ea.	24564
250µL	N	22s	2"/51mm	2	1725	81100	ea.	24567
250µL	RN	22s	2"/51mm	2	1725	81130	ea.	24568
500µL	LTN	22	2"/51mm	2	1750	81217	ea.	24571
500µL	RN	22	2"/51mm	2	1750	81230	ea.	24572
1mL	LTN	22	2"/51mm	2	1001	81317	ea.	24575
1mL	RN	22	2"/51mm	2	1001	81330	ea.	24576
1mL	TLL	*		w/out slots	1001	81320	ea.	24578
2.5mL	TLL	*		w/out slots	1002	81420	ea.	24584
2.5mL	RN	22	2"/51mm	2	1002	81430	ea.	24582
2.5mL	N	22	2"/51mm	2	1002	81417	ea.	24581
5mL	TLL	*		w/out slots	1005	81520	ea.	20178
10mL	TLL	*		w/out slots	1010	81620	ea.	20179
25mL	TLL	*		w/out slots	1025	82520	ea.	24591
50mL	TLL	*		w/out slots	1050	85020	ea.	24717
100mL	TLL	*		w/out slots	1100	86020	ea.	24724

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		Restek	price
					Model	cat.#		
10µL	F	26	50mm	2	10F-GT	002200	ea.	24725
10µL	R	26	50mm	2	10R-GT	002250	ea.	24726
25µL	F	25	50mm	2	25F-GT	003200	ea.	24727
25µL	R	25	50mm	2	25R-GT	003250	ea.	24728
50µL	F	25	50mm	2	50F-GT	004200	ea.	24729
50µL	R	25	50mm	2	50R-GT	004250	ea.	24730
100µL	F	25	50mm	2	100F-GT	005200	ea.	24734
100µL	R	25	50mm	2	100R-GT	005250	ea.	24735
100µL	LL	*			100F-LL-GT	005230	ea.	24737
250µL	F	25	50mm	2	250F-GT	006200	ea.	24738
250µL	R	25	50mm	2	250R-GT	006250	ea.	24739
250µL	LL	*			250F-LL-GT	006230	ea.	24741
500µL	F	25	50mm	2	500F-GT	007200	ea.	24742
500µL	R	25	50mm	2	500R-GT	007250	ea.	24743
1mL	R	23	50mm	2	1MR-GT	008100	ea.	24750
1mL	LL	*			1MDF-LL-GT	008025	ea.	24752
2.5mL	LL	*			2.5MDF-LL-GT	008425	ea.	24755
5mL	LL	*			5MDR-LL-GT	008760	ea.	24757
10mL	LL	*			10MDR-LL-GT	008960	ea.	24759
25mL	LL	*			25MDR-LL-GT	009462	ea.	24760
50mL	LL	*			50MR-LL-GT	009660	ea.	24761
100mL	LL	*			100MR-LL-GT	009760	ea.	24762

*Needles sold separately. See page 382.



24558

24728

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) Teflon® Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip