



21260

24868

### Gas-Tight Syringes for Rheodyne® & Valco® Valves

- Teflon®-tipped plungers.
- Removable needles.
- Replaceable syringe barrels, plungers, and plunger tips.

#### Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	Hamilton cat.#	Restek qty.	Restek cat.#	price
10µL	RN	26s	2"/51mm	3	1701	80065	ea.	21260	
25µL	RN	22s	2"/51mm	3	1702	80265	ea.	21261	
50µL	RN	22s	2"/51mm	3	1705	80965	ea.	21262	
100µL	RN	22s	2"/51mm	3	1710	81065	ea.	21263	
250µL	RN	22s	2"/51mm	3	1725	81165	ea.	21264	

#### SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	SGE cat.#	Restek qty.	Restek cat.#	price
10µL	R	22	2"	3	10R-GT-LC	002313	ea.	24866	
25µL	R	22	2"	3	25R-GT-LC	003312	ea.	24867	
50µL	R	22	2"	3	50R-GT-LC	004312	ea.	24868	
100µL	R	22	2"	3	100R-GT-LC	005312	ea.	24869	
250µL	R	22	2"	3	250R-GT-LC	006312	ea.	24870	
500µL	R	22	2"	3	500R-GT-LC	007312	ea.	24871	

### Replacement Needles for Gas-Tight Syringes for Rheodyne® & Valco® Valves

#### Hamilton



24941

Syringe Volume	Needle Gauge	Needle Length	Point Style	Hamilton cat.#	Restek qty.	Restek cat.#	price
5-100µL	22s	2"	3	7770-01	6-pk.	24941	
250µL-10mL	22	2"	3	7780-04	6-pk.	24945	

#### SGE



24808

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	Restek qty.	Restek cat.#	price
10µL	22	2"	3	037250	5-pk.	24808	
25-500µL	22	2"	3	038250	5-pk.	24809	



## Restek's Learning Network

Sign up for our widely acclaimed seminars today!  
Visit [www.restek.com/seminars](http://www.restek.com/seminars)