



Selection Guide for Q-sep™ dSPE Tubes

Commodity types and examples	AOAC 2007.01	EN 15662	Mini-multiresidue	Additional products
 <p>General purpose</p> <ul style="list-style-type: none"> • Celery • Head lettuce • Cucumber • Melon 	<p>Q-sep Q250 2mL, 100-pk. (cat.# 26124)</p> <p>Q-sep Q350 15mL, 50-pk. (cat.# 26220)</p>	<p>Q-sep Q210 2mL, 100-pk. (cat.# 26215)</p> <p>Q-sep Q370 15mL, 50-pk. (cat.# 26223)</p>	<p>Q-sep Q210 2mL, 100-pk. (cat.# 26215)</p>	
 <p>Fatty or waxy fruits & vegetables</p> <ul style="list-style-type: none"> • Cereals • Avocado • Nuts & seeds • Dairy 	<p>Q-sep Q251 2mL, 100-pk. (cat.# 26125)</p> <p>Q-sep Q351 15mL, 50-pk. (cat.# 26221)</p>		<p>Q-sep Q211 2mL, 100-pk. (cat.# 26216)</p>	<p>Q-sep Q373 15mL, 50-pk. (cat.# 26226)</p>
 <p>Pigmented fruits & vegetables</p> <ul style="list-style-type: none"> • Strawberries • Sweet potatoes • Tomatoes 	<p>Q-sep Q352 15mL, 50-pk. (cat.# 26222)</p>	<p>Q-sep Q212 2mL, 100-pk. (cat.# 26217)</p> <p>Q-sep Q371 15mL, 50-pk. (cat.# 26224)</p>	<p>Q-sep Q212 2mL, 100-pk. (cat.# 26217)</p>	<p>Q-sep Q253 2mL, 100-pk. (cat.# 26123)</p>
 <p>Highly pigmented fruits & vegetables</p> <ul style="list-style-type: none"> • Red peppers • Spinach • Blueberries 	<p>Q-sep Q252 2mL, 100-pk. (cat.# 26219)</p>	<p>Q-sep Q213 2mL, 100-pk. (cat.# 26218)</p> <p>Q-sep Q372 15mL, 50-pk. (cat.# 26225)</p>	<p>Q-sep Q213 2mL, 100-pk. (cat.# 26218)</p>	<p>Q-sep Q374 15mL, 50-pk. (cat.# 26126)</p>
<p>Download free instructions at www.restek.com/quenchers</p>	<p>Instruction sheet# 805-01 002</p>	<p>Instruction sheet# 805-01 001</p>	<p>Instruction sheet# 805-01 001</p>	<p>Generic dSPE 805-01 003</p>



Restek Customer Service

In the U.S.

Call: 800-356-1688 (ext. 3) or 814-353-1300 (ext. 3)


Monday–Friday 8:00 a.m.–6:00 p.m. ET

Fax: 814-353-1309—24-hours a day

Online: www.restek.com—24-hours a day

Outside the U.S.

Contact your Restek representative:
Refer to our list on pages 4-5 or visit our website at www.restek.com



Melissa Decker, Customer Service