



Lid components prior to assembly.



Fully assembled unit shown with Glass Vacuum Chamber (cat.# 25991) and Disk Holders (cat.# 24020)

Resprep® 6-Port Disk Manifold Lid

- Low cost option for disk extraction; fits standard 3 3/4" x 7 1/2" glass vacuum chambers.
- Doubles sample capacity—holds 6 disks, compared to standard 12-port manifolds, which only hold 3.
- Individual vacuum control for each port improves reproducibility.
- Collection plate design secures variety of receiving vessels.
- Inert Teflon® sample guides reduce cross-contamination and carryover.
- Compatible with any standard male luer end disk holder.



Description	qty.	cat.#	price
Resprep 6-Port Disk Manifold Lid*			
Includes: polypropylene lid with 6 ports, 6 nylon valves, 6 Teflon needle guides, 4 black lid legs, collection baseplate, collection plate for volumetric flasks, collection plate for concentrator tubes, 3 white collection plate posts, 12 collection plate retaining clips	kit	25992	
Accessories			
Glass Vacuum Chamber w/gauge & valve for Resprep manifolds, 6 or 12-Port	ea.	25991	
Manifold Lid Replacement Parts			
Valves, 12 or 24-Port	2-pk.	26083	
Valve Attachment	48-pk.	26130	
Needle Attachment	48-pk.	26131	
Sample Guides, 12 or 24-Port	12-pk.	26084	
Gasket, 12-port	2-pk.	24011	

*Vacuum chamber not included.



Resprep® Sample Delivery System

- Compatible with Resprep® 1, 3, 6, and 15 mL SPE cartridges and Diskcover-47 Extraction Disk Holder (cat.# 24020).
- 6 Teflon® transfer lines (1/8" OD x 1/16" ID x 36" long); each contains a different embedded color stripe, for easy sample identification.
- Specified in EPA drinking water methods.
- Tested to pH of 1 to ensure no contaminants leach from system.

Use the Resprep® Sample Delivery System to transfer large volumes of low viscosity samples directly from a bottle to a solid phase extraction cartridge or to a disk on a vacuum manifold system for extraction or cleanup. Each unit consists of 6 transfer lines with a stainless steel weight on one end, and a color-coded screw fitting and polyethylene terephthalate (PET) adapter on the opposite end.

Description	qty.	cat.#	price
Resprep Sample Delivery System	6-pk.	26250	

ChromaBLOGraphy
 Topical and timely insights from top chromatographers.
 Visit us at blog.restek.com