

### Resprep®-C18 and Resprep®-C8 SPE Disks

- Glass fiber disks embedded with C18 or C8 bonded silica.
- Extract semivolatiles organic compounds.
- Deep-pore design reduces clogging and allows faster flow rates.
- Meet requirements for EPA Methods 525.1, 506, 550.1, and 549.1.
- Lower cost than Teflon® disks.

Description	Porosity	qty.	cat.#	price
Resprep-C8	47mm	24-pk.	24048	
Resprep-C18	47mm	20-pk.	24004	
Resprep-C18	90mm	12-pk.	25988	

### Resprep® Oil & Grease SPE Disks

- 47 mm glass fiber disks embedded with specialty bonded silica.
- Meet requirements for EPA Method 1664.\*
- Reduce emulsion formation and amount of solvent required by previous EPA methods.
- No chlorofluorocarbons needed.

Description	qty.	cat.#	price
Resprep Oil & Grease SPE Disks	20-pk.	26022	

\*A sodium sulfate drying tube and a 0.45 µm PTFE syringe filter (cat.# 26145, page 399) also may be used.

### Resprep® SPE Flow Filters

- Designed specifically to improve flow when filtering oil and grease samples.
- Use with Resprep® Diskcover-47 reservoir, or any 47 mm glass sample reservoir.

Description	qty.	cat.#	price
Resprep SPE Flow Filters	20-pk.	26024	

### Resprep® Resin SPE Disks

- 47 mm glass fiber disks embedded with styrene/DVB resin.
- For chlorinated, benzidine-containing, or nitrogen-containing pesticides.
- Meet requirements of EPA Methods 515.2 and 553.

Description	qty.	cat.#	price
Resprep Resin SPE Disks	20-pk.	26023	

## Parts for Diskcover-47 Extraction System

### Diskcover-47 Extraction Disk Holder

- Compatible with most vacuum manifold systems that accept 1/8-inch male luer fittings.
- Sample can be automatically introduced via 1/8-inch Teflon® tubing or from the optional Diskcover-47 reservoir.

Description	qty.	cat.#	price
Diskcover-47 Extraction Disk Holder	ea.	24020	
Diskcover-47 Extraction Disk Holder	6-pk.	24021	
Teflon Tube Luer Adaptors (1/8" OD)	6-pk.	24017	
Teflon Sample Tubing (2ft./0.6m, 1/8" OD)	6-pk.	24025	

### Diskcover-47 Reservoir\*

The Diskcover-47 open-top reservoir allows you to pour up to 125 mL of sample directly onto the filter disk holder. It easily installs on top of the Diskcover-47 extraction disk holder.

Description	qty.	cat.#	price
Diskcover-47 Reservoir	ea.	24029	
Diskcover-47 Reservoir	6-pk.	24030	
Teflon O-Ring**	ea.	24027	
Teflon O-Ring**	6-pk.	24028	

\*Must be used with the Diskcover-47 extraction disk holder.

\*\*Order when using PTFE extraction filters.




Resprep® disks & flow filters extract analytes of interest at high flow rates, and significantly reduce clogging.

also available

For an alternative membrane filtration system, see page 343.



24020



24029