

Accelerated Solvent Extraction (ASE)

Accelerated solvent extraction is a common technique for fast and reliable extraction of organic materials from solid matrices, using EPA SW-846 Method 3545, Pressurized Fluid Extraction (PFE). Restek offers a wide range of replacement parts to keep your extraction system running smoothly.

Replacement Parts for ASE® 200 Systems, Manufacturer's Design



Cell bodies are now serialized for easy sample identification.

Extraction Cell Bodies for ASE® 200 Systems

- Cell bodies are now serialized for easy sample identification.
- Choose original equipment-equivalent stainless steel, or Siltek® deactivation for improved inertness.
- Inner surfaces polished, for easier cleaning.
- Caps include frit, washer, PTFE O-ring, and threaded insert.

Extraction Cell Body	Similar to		Stainless Steel		Siltek Treated		
	Dionex part #	qty.	cat.#	price	qty.	cat.#	price
1mL for ASE 200	054973	ea.	26110		ea.	26111	
5mL for ASE 200	054974	ea.	26112		ea.	26113	
11mL for ASE 200	048820	ea.	26114		ea.	26115	
22mL for ASE 200	048821	ea.	26098		ea.	26099	
33mL for ASE 200	048822	ea.	26116		ea.	26117	

Extraction Cell Replacement Parts for ASE® 200 Systems

Description	Similar to		Stainless Steel		Siltek Treated		
	Dionex part #	qty.	cat.#	price	qty.	cat.#	price
Replacement Extraction Cell End Caps for ASE 200	049450	2-pk.	26096		2-pk.	26097	
Cap Inserts for ASE 200		2-pk.	26166				
Replacement Frits for ASE 200	049453	10-pk.	26100		10-pk.	26101	

Description	Similar to		Stainless Steel		Siltek Treated		
	Dionex part #	qty.	cat.#	price	qty.	cat.#	price
Snap Rings for Caps for ASE 200	049456	10-pk.	26184				
Funnel for ASE 200	056958	ea.	26180				
PTFE O-Rings for ASE 200 & ASE 300 Caps	049457	100-pk.	26187				
Viton O-Rings for ASE 200 & ASE 300 Caps	056325	50-pk.	26188				



26166



26180



26256



26119



26121

26122

PEEK™ Washers for ASE® 200 Systems

- Meet original equipment manufacturer's performance.
- Available in 4 quantities.



Description	Similar to Dionex part #	qty.	cat.#	price
PEEK Washers for ASE 200	049454	12-pk.	25256	
		48-pk.	25257	
		250-pk.	26120	
		1,000-pk.	26229	

20 mm Filters for ASE® 200 Extraction Cells

Description	Similar to Dionex part #	qty.	cat.#	price
Cellulose Filters for ASE 200	049458	100-pk.	26118	
Cellulose Filters for ASE 200	049458	1,000-pk.	26190	
Glass Fiber Filters for ASE 200	047017	100-pk.	26119	

60 mL Sample Collection Vials for ASE® 200 Systems

- Cleaned, assembled, and ready to use.
- Caps and Teflon®-lined septa included.

Description	Color	Similar to Dionex part #	qty.	cat.#	price
60mL Collection Vials for ASE 200	Clear	048784	72-pk.	26121	
60mL Collection Vials for ASE 200	Amber	048781	72-pk.	26122	
Replacement Septa (24mm x 0.125" Teflon-lined silicone)		—	100-pk.	24694	