

Replacement Parts for ASE® 300/100 Systems, Manufacturer's Design

Extraction Cell Bodies for ASE® 300/100 Systems

- Cell bodies are now serialized for easy sample identification.
- Choose original equipment-equivalent stainless steel, or Siltek® deactivation for improved inertness.
- Inner surfaces polished, for easier cleaning.
- Caps include frit, washer, PTFE O-ring, and threaded insert.

Extraction Cell Body	Similar to		Stainless Steel		Siltek Treated		
	Dionex part #	qty.	cat.#	price	qty.	cat.#	price
10mL for ASE 100	059793	ea.	26172		ea.	26173	
34mL for ASE 300/100	056646	ea.	26176		ea.	26177	
66mL for ASE 300/100	056696	ea.	26178		ea.	26179	
100mL for ASE 300/100	056693	ea.	26132		ea.	26133	



Cell bodies are now serialized for easy sample identification.

Extraction Cell Replacement Parts for ASE® 300/100 Systems

Description	Similar to		Stainless Steel		Siltek Treated		
	Dionex part #	qty.	cat.#	price	qty.	cat.#	price
Replacement Extraction Cell End Caps for ASE 300	056921	2-pk.	26170		2-pk.	26171	
Cap Inserts for ASE 300		2-pk.	26167				
Replacement Frits for ASE 300/100		6-pk.	26174		6-pk.	26175	
Description	Similar to		Stainless Steel		Siltek Treated		
	Dionex part #	qty.	cat.#	price	qty.	cat.#	price
Snap Rings for Caps for ASE 300/100	056778	12-pk.	26134				
Funnel for ASE 300/100	056699	ea.	26169				
PTFE O-Rings for ASE 200 & ASE 300 Caps	049457	100-pk.	26187				
Viton O-Rings for ASE 200 & ASE 300 Caps	056325	50-pk.	26188				



PEEK™ Washers for ASE® 300/100 Systems

- Meet original equipment manufacturer's performance.
- Available in 2 quantities.

Description	Similar to Dionex part #	12-pk.	48-pk.
PEEK Washers for ASE 300/100	061687	25393	25394

30 mm Filters for ASE® 300/100 Extraction Cells

Description	Similar to Dionex part #	qty.	cat.#	price
Cellulose Filters for ASE 300/100	056780	100-pk.	26168	
Glass Fiber Filters for ASE 300/100	056781	100-pk.	26189	



250 mL Sample Collection Bottles for ASE® 300/100 Systems

- Cleaned, assembled, and ready to use.
- Caps and Teflon®-lined septa included.

Description	Color	Similar to Dionex part #	qty.	cat.#	price
250mL Collection Bottles for ASE 300/100	Clear Glass	056284	12-pk.	26191	
250mL Collection Bottles for ASE 300/100	Amber Glass	—	12-pk.	26260	
Replacement Septa (24mm x 0.125" Teflon-lined silicone)		—	100-pk.	24694	



also **available**

Diatomaceous Earth and Sodium Sulfate
 See page 398.