

## Flasks


### Erlenmeyer Flasks w/Stopper

ST outer joint.

Description	Volume	Taper Size	qty.	cat.#	price
Erlenmeyer Flask w/Stopper	125mL	24/40	ea.	23345	
Erlenmeyer Flask w/Stopper	250mL	24/40	ea.	23346	
Erlenmeyer Flask w/Stopper	500mL	24/40	ea.	23347	
Erlenmeyer Flask w/Stopper	1,000mL	24/40	ea.	23348	



23345



23350

### Erlenmeyer Flasks w/Stopper, Iodine Style

Full length ST joint, funnel top and hollow stopper with closed bottom.


Description	Volume	Taper Size	qty.	cat.#	price
Erlenmeyer Flask w/Stopper, Iodine Style	125mL	24/40	kit	23349	
Erlenmeyer Flask w/Stopper, Iodine Style	250mL	24/40	kit	23350	
Erlenmeyer Flask w/Stopper, Iodine Style	500mL	24/40	kit	23351	\$54

### Volumetric Flasks, Class A, ST Stopper

Calibrated to contain, with marking spot and ST stopper. Precision grade and “unserialized” flasks which meet the requirements for type 1, Class A, Federal Specification NNN-F-289 and ASTM Specification E 288.

Description	Volume	Taper Size	qty.	cat.#	price
Volumetric Flask w/Stopper	1mL*	ST# 9	ea.	23359	
Volumetric Flask w/Stopper	5mL	ST# 9	ea.	23360	
Volumetric Flask w/Stopper	10mL	ST# 9	ea.	23362	
Volumetric Flask w/Stopper	25mL	ST# 9	ea.	23363	
Volumetric Flask w/Stopper	50mL	ST# 9	ea.	23364	
Volumetric Flask w/Stopper	100mL	ST# 13	ea.	23365	
Volumetric Flask w/Stopper	200mL	ST# 16	ea.	23366	
Volumetric Flask w/Stopper	250mL	ST# 16	ea.	23367	
Volumetric Flask w/Stopper	500mL	ST# 19	ea.	23368	
Volumetric Flask w/Stopper	1000mL	ST# 22	ea.	23369	

\*The 1mL is test tube shaped and not guaranteed to the Class A accuracy.




23365

### Volumetric Flasks, Class A, ST Stopper, Amber

Calibrated to contain, with marking spot and ST stopper. Amber stain process provides protection to light-sensitive liquids. Meet the requirements for type 1, Class A, Federal Specification NNN-F-289 and ASTM Specification E 288.

Description	Volume	Taper Size	qty.	cat.#	price
Volumetric Flask w/Stopper, Amber	10mL	ST# 9	ea.	23352	
Volumetric Flask w/Stopper, Amber	25mL	ST# 9	ea.	23353	
Volumetric Flask w/Stopper, Amber	50mL	ST# 9	ea.	23354	
Volumetric Flask w/Stopper, Amber	100mL	ST# 13	ea.	23355	
Volumetric Flask w/Stopper, Amber	250mL	ST# 16	ea.	23356	
Volumetric Flask w/Stopper, Amber	500mL	ST# 19	ea.	23357	
Volumetric Flask w/Stopper, Amber	1000mL	ST# 22	ea.	23358	



23355