



23374

Flasks *cont'd*

Micro Volumetric Flask, Class A, ST Stopper

The 1 through 5 mL sizes conform to the recommendations of the Committee on Microchemical Apparatus, Analytical Division of the A.C.S., in their 1956 report. Sizes 10 through 50 mL have necks of special design, to allow pipetting to the very bottom of the flask. The wide base gives much greater support to the flask.

Description	Volume	Taper Size	qty.	cat.#	price
Micro Volumetric Flask w/Stopper	1mL	ST# 8	kit	23370	
Micro Volumetric Flask w/Stopper	2mL	ST# 8	kit	23372	
Micro Volumetric Flask w/Stopper	5mL	ST# 8	kit	23373	
Micro Volumetric Flask w/Stopper	10mL	ST# 13	kit	23374	
Micro Volumetric Flask w/Stopper	25mL	ST# 13	kit	23375	
Micro Volumetric Flask w/Stopper	50mL	ST# 13	kit	23376	



23377



Squeeze-O-Matic Burets

Burets

Auto-Zero Burets with Teflon® Stopcock and Reservoir, Class A

Automatic, precision grade white graduations. Supplied with bottle and bulb.

Description	Volume	qty.	cat.#	price
Auto-Zero Buret w/1,000mL Reservoir	25mL	kit	23377	
Auto-Zero Buret w/2,000mL Reservoir	50mL	kit	23378	

Squeeze-O-Matic (S-O-M) Buret with Poly Bottle and Teflon® Plug

Description	Volume	Color	qty.	cat.#	price
S-O-M Buret, PTFE Plug	10mL	white	kit	23379	
S-O-M Buret, PTFE Plug	25mL	white	kit	23380	



23381

Funnels

Separatory Funnel, Squibb Type with Teflon® Stopcock

With solid ST stopper.

Description	Volume	Stopcock	qty.	cat.#	price
Separatory Funnel, Squibb Type	30mL	2 PTFE	kit	23381	
Separatory Funnel, Squibb Type	60mL	2 PTFE	kit	23382	
Separatory Funnel, Squibb Type	125mL	2 PTFE	kit	23383	
Separatory Funnel, Squibb Type	250mL	4 PTFE	kit	23384	
Separatory Funnel, Squibb Type	500mL	4 PTFE	kit	23385	
Separatory Funnel, Squibb Type	1L	4 PTFE	kit	23386	
Separatory Funnel, Squibb Type	2L	6 PTFE	kit	23387	