

Glass Impingers

Midget Glass Impingers

Use with a sampling pump to trap air contaminants into liquid collection media, as specified in OSHA and NIOSH industrial hygiene methods. Both dispersion and fritted nozzles are available as bubblers.

Description	Volume	Taper Size	qty.	cat.#	price
Midget Glass Impinger w/Fritted Tube	30mL	24/40	kit	23388	
Midget Glass Impinger w/Dispersion Tube	30mL	24/40	kit	23389	



Pipets

Automatic, Repeating Pipets

Delivers identical amounts of liquid in rapid succession, by tilting dispenser to the rear, allowing reagent to fill the calibrated section of pipet from Erlenmeyer flask. The head and flask are held together by ST joints and the included poly joint clamp. Accuracy within 2%.

Description	Volume	Taper Size	qty.	cat.#	price
Repeating Pipet w/250mL Reservoir	1mL	24/40	kit	23390	
Repeating Pipet w/250mL Reservoir	5mL	24/40	kit	23391	
Repeating Pipet w/500mL Reservoir	10mL	24/40	kit	23392	
Repeating Pipet w/500mL Reservoir	15mL	24/40	kit	23393	
Repeating Pipet w/500mL Reservoir	20mL	24/40	kit	23394	
Repeating Pipet w/1L Reservoir	25mL	24/40	kit	23395	



23392

Stopper

ST Hollow Stopper

Pennyhead, full length glass apparatus stopper with closed bottom.

Description	Taper Size	qty.	cat.#	price
ST Hollow Stopper	24/40	ea.	23338	



23338



Catch the Buzz

Sign up for Restek's e-newsletter, *The Buzz*
www.restek.com/buzz