

Buy only the parts you need!

Replacement Orifices

Use these orifices with a Veriflo® 423XL flow controller to change the flow range for alternative sampling times.

Flow	Orifice size	Siltek Treated		Stainless Steel	
		cat.#	price	cat.#	price
0.5–2 mL/min.	0.0008"	24219		24218	
2–4 mL/min.	0.0012"	24233		24245	
4–8 mL/min.	0.0016"	24234		24246	
8–15 mL/min.	0.0020"	24235		24247	
15–30 mL/min.	0.0030"	24236		24248	
30–80 mL/min.	0.0060"	24237		24249	
80–340 mL/min.	0.0090"	22099		22098	



Critical orifice



24249

2 µm Frit Filters

For use in critical orifice fitting. Includes washers.

Description	qty.	cat.#	price
Replacement Frit Filter, Siltek Treated	3-pk.	24171	
Replacement Frit Filter, Stainless Steel	3-pk.	24170	



24171


24170

Veriflo® Flow Controllers

Veriflo® 423XL flow controllers are offered in a Siltek® and stainless steel version. The flow device is available with or without a critical orifice. (Vacuum gauge sold separately.)

The critical orifice in a Veriflo® flow controller is interchangeable. Order orifices for alternate sampling times, or replacement orifices, separately.

Flow	Orifice size	Siltek Treated		Stainless Steel	
		cat.#	price	cat.#	price
0.5–2 mL/min.	0.0008"	24232		24229	
2–4 mL/min.	0.0012"	24255		24260	
4–8 mL/min.	0.0016"	24256		24261	
8–15 mL/min.	0.0020"	24257		24262	
15–30 mL/min.	0.0030"	24258		24263	
30–80 mL/min.	0.0060"	24259		24264	
80–340 mL/min.	0.0090"	22103		22102	
	without orifice	24238		24239	



Flow controller



24262

7µm In-Line Filter

This 316 stainless steel filter is designed to collect particles larger than 7 microns. We offer Siltek® and stainless steel versions (1/4" compression fitting on both ends).

Description	qty.	cat.#	price
7µm In-Line Filter, Siltek Treated	ea.	24265	
7µm In-Line Filter, Stainless Steel	ea.	24266	

Note: frit is not replaceable.



24266

Sample Inlets

- 1/4" stainless steel compression fitting on each end.
- One end connects to flow controller or canister; nut on other end serves as rain guard.
- Includes nuts and ferrules.
- Two different lengths for use with large canisters and miniature canisters.

Description	qty.	Siltek Treated		Stainless Steel	
		cat.#	price	cat.#	price
Sample Inlet, 6" Length	ea.	26210		26209	
Sample Inlet, 1.5" Length	ea.	26212		26211	



26209

26211