

Tedlar® Sampling Bags

- Find the bags you need—we offer sizes from 0.5 liters to 100 liters.
- Unique all-in-one septum and valve fitting make these lightweight and easy to use.
- Polypropylene or stainless steel valve.
- Both valves conveniently connect to 3/16" ID Teflon® tubing.

Description	qty.	Polypropylene Valve		Stainless Steel Valve	
		cat.#	price	cat.#	price
0.5L 6" x 6"	10-pk.	22049		22038	
1L 7" x 7"	10-pk.	22050		22039	
3L 9.5" x 10"	10-pk.	22051		22040	
5L 12" x 12.5"	10-pk.	22052		22041	
10L 11.75" x 22"	10-pk.	22053		22042	
12L 13" x 24"	10-pk.	22054		22043	
25L 17.5" x 24"	5-pk.	22055		22044	
40L 24" x 24.25"	5-pk.	22056		22045	
80L 28.25" x 30.5"	5-pk.	22057		22046	
100L 28" x 36"	3-pk.	22058		22047	

Description	qty.	cat.#	price
Teflon Faced Silicone Replacement Septum, 4mm diameter	10-pk.	22104	



tech tip

Use septum puller (cat.# 20117) to replace septum in sampling bag valve.



Description	qty.	cat.#	price
Septum Puller	ea.	20117	\$16.25

also available

Multi-layer foil bags

Visit www.restek.com

Vacuum Bag Sampler

- Fast bag sampling without sample passing through pump.
- Bag capacity up to 10 L.

The Model 1062 Vacuum Bag Sampler provides fast sampling with zero cross-contamination. A vacuum created in the box draws air into the sampling bag without drawing it through the vacuum pump first, as is the case with standard air sampling pumps, thereby preventing contamination of the sample. This bag sampler can fill a 10 L bag in two minutes with an automatic shut-off switch, which stops the sample bag from overfilling. The filling rate is adjusted with a vent rotometer valve. An external battery recharging port enables continuous operation with battery charger. In addition, the quick exhaust valve allows for fast removal of the sampling bag. The sampler comes with a universal power adapter/charger, battery, instruction manual, and 1-year limited warranty.

Specifications:

Sampling Bag:	1 bag up to 10L size
Running Time:	8 hours
Flow Rate (Fill Rate):	1-5L/min.
Power Requirements:	12V battery, 4.5 amp
Charge Time:	9 hours
Dimensions:	9" x 14.6" x 21.7"
Weight:	17 lbs

Description	qty.	cat.#	price
Vacuum Bag Sampler Model 1062 (includes: power adapter, battery, manual)	ea.	22118	
Replacement Battery for Vacuum Bag Sampler Model 1062	ea.	22119	
Universal Battery Charger for Vacuum Bag Sampler Model 1062 (115/230 VAC)	ea.	22120	



Features:

- Observation window on case lid
- Sample inlet accepts 1/4" OD tubing
- Case designed for rugged outdoor use
- CE certified

ChromaBLOGraphy
 Topical and timely insights from top chromatographers.
 Visit us at blog.restek.com