

Gas Regulators for Transportable Cylinders

For this cylinder:

DCG Partnership Cylinders:
Size: 7.6 x 24 cm
CGA-170/110 connection.
US DOT Specs: DOT-4B-240ET
Please note: This cylinder is not approved for use in Canada.



Use this regulator:

Mini-Regulator for natural gas and refinery gas standards

- 0–300 psig inlet pressure range.
- 0–15 psig outlet pressure range.
- Supplied with 0–15 psig outlet pressure gauge, brass CGA 170 nut and nipple.



22032

Description	qty.	cat.#	price
Mini-Regulator	ea.	22032	

For these cylinders:

Spectra (Linde) 104L:
 Aluminum construction
 Size: 8 x 24 cm
 Volume/Pressure:
 104 liters of gas
 @ 1,800 psi
 CGA-180 outlet fitting.
 Weight: 1.5 lbs/0.7 kg



Scotty® (Air Liquide) 110L (Pi-marked Cylinders for EU Regulations):
 Aluminum construction
 Size: 8.3 x 29.5 cm
 Volume/Pressure:
 110 liters of gas @ 1,800 psi
 CGA-180 outlet fitting.
 Weight: 2.2 lbs/1 kg
 DOT Specifications: 3AL2216



Use these regulators:

Spectra Gas 7621 High-Purity VOC Regulator

- Single-stage, stainless steel.
- Two pressure gauges and CGA-180 fitting.
- 3,000 psig maximum inlet pressure.
- Stainless steel diaphragm and Kel-F® seat.
- 1/8-inch tube compression outlet.
- Low internal volume: 3.03 cc.
- Accurate pressure control even at low flow rates.
- Individually tested for leaks and impurities.



21572

Description	qty.	cat.#	price
0–30psig outlet pressure gauge	ea.	21572	
0–100psig outlet pressure gauge	ea.	21572-R100	

See next page for a syringe adapter kit.

Continued on next page.