

Gas Sampling Valves and Sample Loops (Sulfinert® Treated)

- Ideal for samples containing low concentrations of sulfur or other active compounds.
- Sample loop sizes from 5 µL to 5 cc.

Sulfinert® treatment eliminates active sites in the valve or loop, for better recovery of active compounds.

Gas Sampling Valves

(1/16" Fittings, 0.40 mm Port Diameter; "W Type" Valve)

Description	qty.	cat.#	price
Sulfinert Gas Sampling Valve; 4-Port	ea.	20584	
Sulfinert Gas Sampling Valve; 6-Port	ea.	20585	
Sulfinert Gas Sampling Valve; 10-Port	ea.	20586	



20585

Replacement Rotors (Not Coated)

Description	qty.	cat.#	price
Replacement Rotor (not coated) for 4-Port Sulfinert Gas Sampling Valve	ea.	20587	
Replacement Rotor (not coated) for 6-Port Sulfinert Gas Sampling Valve	ea.	20588	
Replacement Rotor (not coated) for 10-Port Sulfinert Gas Sampling Valve	ea.	20589	

Gas Sample Loops

(1/16" fittings, for "W Type" valves)

Description	Size	qty.	cat.#	price
Sample Loops, Sulfinert Treated	5µL	ea.	22840	
Sample Loops, Sulfinert Treated	10µL	ea.	22841	
Sample Loops, Sulfinert Treated	20µL	ea.	22842	
Sample Loops, Sulfinert Treated	25µL	ea.	22843	
Sample Loops, Sulfinert Treated	50µL	ea.	22844	
Sample Loops, Sulfinert Treated	100µL	ea.	22845	
Sample Loops, Sulfinert Treated	250µL	ea.	22846	
Sample Loops, Sulfinert Treated	500µL	ea.	22847	
Sample Loops, Sulfinert Treated	1mL	ea.	22848	
Sample Loops, Sulfinert Treated	2mL	ea.	22849	
Sample Loops, Sulfinert Treated	5mL	ea.	22850	



Jumbo Syringe

Clear acrylic syringes, ideal for holding and dispensing large volumes of gas. An adjustable plunger on the O-ring ensures that the syringe is gas-tight over a long period of time. The central port is supplied with a luer-lock fitting; the secondary port is supplied with a septum nut. This enables access to the gas sample for adding standards or removing a subsample. The plunger stem is detachable, making sample storage easy.



21276

Volume	Model	SGE cat.#	qty.	Restek cat.#	price
500mL	500MAR-LL-GT	009910	ea.	21275	
1000mL	1000MAR-LL-GT	009920	ea.	21276	
2000mL	2000MAR-LL-GT	009930	ea.	21277	

Syringe O-Rings

Syringe Volume	SGE cat.#	qty.	Restek cat.#	price
500mL	032527	ea.	21278	
1000mL	032532	ea.	21279	



21279

21278



For more information on all air monitoring canisters and products, visit www.restek.com/air

