

8000 Series Methods

Method 8021 (Volatile Organics)

502.2 Internal Standard Mix #2 (2 components)

2-bromo-1-chloropropane fluorobenzene
2,000µg/mL each in P&T methanol, 1mL/ampul
cat. # 30041 (ea.)

8021 Surrogate Mix (2 components)

2-bromochlorobenzene 1,4-dichlorobutane
1,500µg/mL each in P&T methanol, 1mL/ampul
cat. # 30086 (ea.)

8021/502.2 Surrogate Mix #1 (3 components)

1-bromo-2-chloroethane fluorobenzene
1-chloro-3-fluorobenzene
2,000µg/mL each in P&T methanol, 1mL/ampul
cat. # 30463 (ea.)

8021/502.2 Surrogate Mix #2 (4 components)

1-bromo-2-chloroethane 1-chloro-3-fluorobenzene
4-bromochlorobenzene fluorobenzene
2,000µg/mL each in P&T methanol, 1mL/ampul
cat. # 30464 (ea.)

Method 8040 (Phenols)

8040 Surrogate Mix (2 components)

2-fluorophenol 2,4,6-tribromophenol
2,000µg/mL each in isopropanol, 1mL/ampul
cat. # 31090 (ea.)

8040 Phenols Mix #1 (9 components)

4-chloro-3-methylphenol 4-nitrophenol
2,4-dichlorophenol pentachlorophenol
2-methyl-4,6-dinitrophenol phenol
3-methylphenol 2,4,6-trichlorophenol
2-nitrophenol
2,000µg/mL each in isopropanol, 1mL/ampul
cat. # 31088 (ea.)

8040 Phenols Mix #2 (9 components)

sec-butyl-4,6-dinitrophenol (dinoseb) 2-methylphenol
2-chlorophenol 4-methylphenol
2,6-dichlorophenol 2,3,4,6-tetrachlorophenol
2,4-dimethylphenol 2,4,5-trichlorophenol
2,4-dinitrophenol
2,000µg/mL each in isopropanol, 1mL/ampul
cat. # 31089 (ea.) \$33

Method 8061A (Phthalate Esters)

8061A Matrix Spike Solution (2 components)

benzyl butyl phthalate bis(2-ethylhexyl)phthalate
2,000µg/mL each in acetone, 1mL/ampul
cat. # 31846 (ea.)

Benzyl Benzoate (Internal Standard)

5,000µg/mL in hexane, 1mL/ampul
cat. # 31847 (ea.)

Method 8061A (Phthalate Esters) *cont'd*

8061A Surrogate Standard (3 components)

dibenzyl phthalate diphenyl phthalate
diphenyl isophthalate
500µg/mL each in acetone, 1mL/ampul
cat. # 31848 (ea.)

EPA 8061A Phthalate Esters Mixture (15 components)

benzyl butyl phthalate di-*n*-hexyl phthalate
bis(2-*n*-butoxyethyl)phthalate dimethylphthalate
bis(2-ethoxyethyl)phthalate di-nonyl phthalate
bis(2-ethylhexyl)phthalate di-*n*-octyl phthalate
bis(2-methoxyethyl)phthalate dipentylphthalate
bis(4-methyl-2-pentyl)phthalate phthalic acid dicyclohexyl ester
di-*n*-butylphthalate phthalic acid diisobutyl ester
diethylphthalate
1,000µg/mL each in hexane:acetone (80:20), 1mL/ampul
cat. # 33227 (ea.)

EPA 8061A Hexyl 2-Ethylhexyl Phthalate Standard

hexyl 2-ethylhexyl phthalate
1,000µg/mL in hexane:acetone (80:20), 1mL/ampul
cat. # 33228 (ea.)

EPA 8061A Phthalate Esters Kit

Contains 1mL each of these mixtures.
33227: EPA 8061A Phthalate Esters Mixture
33228: EPA 8061A Hexyl 2-Ethylhexyl Phthalate Standard
cat. # 33229 (kit)



Quantity discounts not available.

Method 8080, 8081 (Chlorinated Pesticides)

Organochlorine Pesticide Mix AB #1 (20 components)

aldrin dieldrin
α-BHC endosulfan I
β-BHC endosulfan II
δ-BHC endosulfan sulfate
γ-BHC (lindane) endrin
α-chlordane endrin aldehyde
γ-chlordane endrin ketone
4,4'-DDD heptachlor
4,4'-DDE heptachlor epoxide (isomer B)
4,4'-DDT methoxychlor
200µg/mL each in hexane:toluene (1:1), 1mL/ampul
cat. # 32291 (ea.)

Get More!

Environmental
Solutions Online

www.restek.com/enviro

