

CLP Pesticides Mixtures, QA Mixes

Pesticide Matrix Spike Mix (6 components)

aldrin	25µg/mL	dieldrin	50
γ-BHC (lindane)	25	endrin	50
4,4'-DDT	50	heptachlor	25
In acetone, 1mL/ampul			
cat. # 32018 (ea.)			

Pesticide Matrix Spike Mix (2/88 SOW) (6 components)

aldrin	200µg/mL	dieldrin	500
γ-BHC (lindane)	200	endrin	500
4,4'-DDT	500	heptachlor	200
In methanol, 1mL/ampul			
cat. # 32031 (ea.)			

Pesticide Evaluation Mix (2/88 SOW) (4 components)

aldrin		4,4'-DDT	
dibutylchloroendate (SS)		endrin	
100µg/mL each in hexane, 1mL/ampul			
cat. # 32032 (ea.)			

See complete listing of PCBs, pages 495-496.

Low Concentration Pesticides Mixtures, QA Mixes

L/C Pesticide Lab Control Sample (7 components)

γ-BHC (lindane)	10µg/mL	endosulfan sulfate	20
γ-chlordane	10	endrin	20
4,4'-DDE	20	heptachlor epoxide (isomer B)	10
dieldrin	20		
In acetone, 1mL/ampul			
cat. # 32040 (ea.)			

CLP Pesticides Mixtures, Calibration Mixes

Pesticide Standard Mix A (2/88 SOW) (10 components)

aldrin	10µg/mL	endosulfan II	20
γ-BHC (lindane)	5	endrin aldehyde	25
4,4'-DDT	20	heptachlor	10
dieldrin	10	heptachlor epoxide (isomer B)	10
endosulfan I	10	methoxychlor	100
In hexane, 1mL/ampul			
cat. # 32033 (ea.)			

Pesticide Standard Mix B (2/88 SOW) (11 components)

aldrin	10µg/mL	4,4'-DDD	20
α-BHC	5	4,4'-DDE	10
β-BHC	10	endosulfan sulfate	20
δ-BHC	10	endrin	10
α-chlordane	10	endrin ketone	20
γ-chlordane	10		
In hexane, 1mL/ampul			
cat. # 32034 (ea.)			

Pesticide Standard Mix A (9 components)

α-BHC	5µg/mL	endosulfan I	5
γ-BHC (lindane)	5	endrin	10
4,4'-DDD	10	heptachlor	5
4,4'-DDT	10	methoxychlor	50
dieldrin	10		
In hexane:toluene (90:10), 1mL/ampul			
cat. # 32297 (ea.)			

CLP Pesticides Mixtures, Calibration Mixes *cont'd*

Pesticide Standard Mix B (11 components)

aldrin	5µg/mL	endosulfan II	10
β-BHC	5	endosulfan sulfate	10
δ-BHC	5	endrin aldehyde	10
α-chlordane	5	endrin ketone	10
γ-chlordane	5	heptachlor epoxide (isomer B)	5
4,4'-DDE	10		
In hexane:toluene (90:10), 1mL/ampul			
cat. # 32298 (ea.)			

Organochlorine Pesticide Mix AB #1 (20 components)

aldrin		dieldrin	
α-BHC		endosulfan I	
β-BHC		endosulfan II	
δ-BHC		endosulfan sulfate	
γ-BHC (lindane)		endrin	
α-chlordane		endrin aldehyde	
γ-chlordane		endrin ketone	
4,4'-DDD		heptachlor	
4,4'-DDE		heptachlor epoxide (isomer B)	
4,4'-DDT		methoxychlor	
200µg/mL each in hexane:toluene (1:1), 1mL/ampul			
cat. # 32291 (ea.)			

Organochlorine Pesticide Mix AB #2 (20 components)

aldrin	8µg/mL	dieldrin	16
α-BHC	8	endosulfan I	8
β-BHC	8	endosulfan II	16
δ-BHC	8	endosulfan sulfate	16
γ-BHC (lindane)	8	endrin	16
α-chlordane	8	endrin aldehyde	16
γ-chlordane	8	endrin ketone	16
4,4'-DDD	16	heptachlor	8
4,4'-DDE	16	heptachlor epoxide (isomer B)	8
4,4'-DDT	16	methoxychlor	80
In hexane:toluene (1:1), 1mL/ampul			
cat. # 32292 (ea.)			

Organochlorine Pesticide Mix AB #3 (20 components)

aldrin		dieldrin	
α-BHC		endosulfan I	
β-BHC		endosulfan II	
δ-BHC		endosulfan sulfate	
γ-BHC (lindane)		endrin	
α-chlordane		endrin aldehyde	
γ-chlordane		endrin ketone	
4,4'-DDD		heptachlor	
4,4'-DDE		heptachlor epoxide (isomer B)	
4,4'-DDT		methoxychlor	
2,000µg/mL each in hexane:toluene (1:1), 1mL/ampul			
cat. # 32415 (ea.)			



also **available**

Use our Rtx®-CLPesticides/Rtx®-CLPesticides2 capillary column pair to analyze all twenty organochlorine pesticides and their surrogates simultaneously. See **pages 88-89**.

See **page 496** for our complete listing of PCBs and congeners.