

Petroleum Standards

Sulfur in Isooctane Calibration Kits and Check Standards *cont'd*

Cal Kit SISO 2.5 - 50 ppm

Blank
 2.50 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 5.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 10.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 15.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 20.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 25.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 50.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 Calibration kit for total sulfur by weight from di-*n*-butylsulfide in isooctane
 2.5–50 ppm range. Set of eight 1mL pre-scored ampuls.
 cat. # 33037 (kit)

Quantity discounts not available.

Check Standard SISO 30

30 ppm total sulfur by weight from di-*n*-butylsulfide in isooctane.
 Set of five 1mL pre-scored ampuls.
 cat. # 33038 (ea.)

Quantity discounts not available.

Cal Kit SISO 50 - 1000 ppm

Blank
 50 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 75 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 100 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 250 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 500 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 1000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 Calibration kit for total sulfur by weight from di-*n*-butylsulfide in isooctane
 50–1000 ppm range. Set of seven 1mL pre-scored ampuls.
 cat. # 33039 (kit)

Quantity discounts not available.

Check Standard SISO 300

300 ppm total sulfur by weight from di-*n*-butylsulfide in isooctane.
 Set of five 1mL pre-scored ampuls.
 cat. # 33040 (ea.)

Quantity discounts not available.

Cal Kit SISO 1000 - 6000

Blank
 1000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 1500 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 2000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 4000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 6000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
 Calibration kit for total sulfur by weight from di-*n*-butylsulfide in isooctane
 1000–6000 ppm range. Set of six 1mL pre-scored ampuls.
 cat. # 33041 (kit)

Quantity discounts not available.

Check Standard SISO 3000

3000 pm total sulfur by weight from di-*n*-butylsulfide in isooctane.
 Set of five 1mL pre-scored ampuls.
 cat. # 33042 (ea.)

Quantity discounts not available.

also available

Custom total sulfur & total nitrogen in isooctane check standards also available.

Call our Analytical Reference Department at 800-356-1688 or 814-353-1300, or your Restek representative for details.

Total Sulfur & Total Nitrogen in Isooctane Calibration Kits

ASTM Methods D3120, D4045, D4629, D5453, D5762, D6069, D6920

Cal Kit SNISO 0.125 - 5.0

Blank
 0.125 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 0.25 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 0.50 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 1.00 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 2.50 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 5.00 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane

Set of seven 1mL pre-scored amber ampuls.

Quantity discounts not available.

Cal Kit SNISO 5.0 - 50.0

Blank
 5.00 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 10.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 25.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 50.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane

Set of five 1mL pre-scored amber ampuls.

cat. # 33058 (kit)

Quantity discounts not available.

Cal Kit SNISO 50.0 - 1000

Blank
 50.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 75.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 100 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 250 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 500 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
 1000 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane

Set of seven 1mL pre-scored amber ampuls.

cat. # 33059 (kit)

Quantity discounts not available.



Documentation Search

Find what you need fast at www.restek.com/documents

- Material safety data sheets
- Certificates of analysis
- Datapacks (by catalog number and/or lot number)