

# Easily detect 1 µg/g melamine with our complete kit for GC/MS

Full GC/MS method kit available with instructions!

Keep up with advances in melamine analysis at [www.restek.com/melamine](http://www.restek.com/melamine)



## Melamine Analysis Kit

Kit includes:

**Column:**

Rxi-5Sil MS w/5 meter Integra-Guard

**Standards:**

- |   |              |
|---|--------------|
| 33247: 1mL Melamine Stock Standard          | (1,000µg/mL) |
| 33248: 1mL Cyanuric Acid Stock Standard     | (1,000µg/mL) |
| 33249: 1mL Ammelide Stock Standard          | (1,000µg/mL) |
| 33250: 1mL Ammeline Stock Standard          | (1,000µg/mL) |
| 33251: 1mL Benzoguanamine Internal Standard | (1,000µg/mL) |
| 33253: 1mL Melamine Mix Standard            | (1,000µg/mL) |

**Derivatization Reagent:**

35607: BSTFA w/1% TMCS, 25g vial

**Accessories:**

50mL centrifuge tubes, 5-pk.  
13mm, 0.45µm nylon syringe filters, 5-pk.

**Easy-to-follow instructions** with procedural check lists to assist with laboratory documentation.

cat. # 33254 (kit)

Quantity discounts not available.

## Melamine and Related Analogs Stock Standard (4 components)

- |  |               |
|--|---------------|
| ammelide   | cyanuric acid |
| ammeline   | melamine      |
| 1,000µg/mL each in diethylamine:water (20:80), 1mL/ampul |               |
| cat. # 33253 (ea.)                                       |               |



## Melamine Stock Standard

- melamine  
1,000µg/mL in diethylamine:water (20:80), 1mL/ampul  
cat. # 33247 (ea.)

## Cyanuric Acid Stock Standard

- cyanuric acid  
1,000µg/mL in diethylamine:water (20:80), 1mL/ampul  
cat. # 33248 (ea.)

## Ammelide Stock Standard

- ammelide  
1,000µg/mL in diethylamine:water (20:80), 1mL/ampul  
cat. # 33249 (ea.)

## Ammeline Stock Standard

- ammeline  
1,000µg/mL in diethylamine:water (20:80), 1mL/ampul  
cat. # 33250 (ea.)

## Benzoguanamine Internal Standard

- benzoguanamine  
1,000µg/mL in pyridine, 1mL/ampul  
cat. # 33251 (ea.)  
1,000µg/mL in pyridine, 5mL/ampul  
cat. # 33252 (ea.)