

GC Chromatograms by Column Phase (cont'd)

Rxi®-55il MS
 amphetamines675
 benzodiazepines682
 cannabinoids674
 cocaine & metabolites673
 dietary supplements684, 685
 dioxins & furans95
 opiates, extracted from blood674
 organochlorine pesticides639
 polycyclic aromatic hydrocarbons97, 555, 556, 636
 semivolatile organics43, 87, 561 - 566
 sympathomimetic amines675

Rxi®-17
 fragrance allergens61

Rxi®-175il MS
 benzofluoranthenes558
 polycyclic aromatic hydrocarbons45, 73, 98, 557, 558, 637, 638

Rxi®-6245il MS
 basic compounds46
 potential genotoxic impurities (PGIs)689, 690
 residual solvents83, 694, 695
 volatile organics103, 547

Stabilwax®
 acrylamide639
 amines, alcohols, chlorides706
 essential oils626, 627
 FAMES59, 621, 646
 flavor volatiles628
 glycols602
 residual solvents84, 697, 698

Stabilwax®-DA
 cresylic acids707
 fatty acids (free)615, 616
 flavor volatiles67, 629, 631

Stabilwax®-DB
 amines66, 704 - 706

Packed/Micropacked Columns

CarboBlack B
 blood alcohol667
 flavor volatiles630

D3606 Application Column
 gasoline128, 656

HayeSep® Q
 solvents702

5% Krytox®/CarboBlack B
 chlorofluorocarbons706

Molesieve 5A
 permanent gases647

Molesieve 13X
 permanent gases647

Refinery Gas Column Set
 hydrocarbon gases129

Res-Sil® C Packing
 OPN129
 n-Octane129

Rt®-1200/Bentone® 34
 aromatics (BTEX)655

Rt®-XLSulfur
 sulfur compounds131, 664
 sulfurs in beverages632

Rtx®-1 on Chromosorb® WAW
 hydrocarbons (motor oil, aviation gas)656

Rtx®-1 SimDist 2887
 simulated distillation661, 662

ShinCarbon ST
 chlorofluorocarbons649
 hydrocarbon gases647, 652
 permanent gases130, 647

10% TCEP
 aromatics, aliphatics654

20% TCEP
 hydrocarbons (motor oil, aviation gas)656

PLOT Columns

Rt®-Alumina BOND/CFC
 chlorinated fluorocarbons74

Rt®-Alumina BOND/KCI
 hydrocarbon gases650, 651, 653

Rt®-Alumina BOND/Na:SO:
 hydrocarbon gases649, 652, 653
 natural gas651

Rt®-Msieve 5A PLOT
 permanent gases109, 648

Rt®-Q-BOND
 solvents111, 701

Rt®-QS-BOND
 hydrocarbon gases652, 654
 ketones703
 solvents702

Rt®-S-BOND
 solvents701

Rt®-U-BOND
 aldehydes703

GC & HPLC Chromatograms
Compound Index
 starting on page 752

EPA Methods Chart

Method #	Method Title	Page #
502.2	Volatile Organic Compounds	595
504.1	1,2-dibromoethane & 1,2-dibromo-3-chloropropane	552
505	Organohalide Pesticides & PCBs	571
506	Semivolatiles Organics	554
508.1	Chlorinated Pesticides	572
524.2	Volatile Organics	551
525.2	Semivolatiles Organics	568
526	Semivolatiles Organics (screening)	560
528	Phenols	560
551.1	Chlorinated Disinfection Byproducts	589, 590
552.2	Haloacetic Acids	603, 604
604	Phenols	553
610	Polycyclic Aromatic Hydrocarbons	555-559
615	Chlorophenoxyacid Herbicides	580
619	Nitrogen/Phosphorus Herbicides & Pesticides	578, 581
624	Purgeable Halocarbons	467
1613	Dioxins	591
1614	Brominated Flame Retardants	588
8021	Volatile Organics	101
8081	Chlorinated Pesticides	53, 569, 570
8082	PCBs	584-587
8095	Explosives	594
8100	Polycyclic Aromatic Hydrocarbons	557
8140/8141/8141A	Organophosphorus Pesticides	573-575
8151A	Chlorophenoxyacid Herbicides	579
8260	Volatile Organics	547-550
8270	Semivolatiles Organics	561-567
8270 Appendix IX	Semivolatiles Organics	566
	TO-14	606, 608
	TO-15	608, 609

Looking for more?

Visit our online technical library for chromatograms based on ASTM methods, USP procedures, and other applications.

www.restek.com/chromatograms

