




**Permanent Gases**  
**ShinCarbon ST**  
(micropacked)




**Permanent Gases, C1+C2 hydrocarbons**  
**ShinCarbon ST**  
(micropacked)



**Permanent Gases**  
**Molecular Sieve 5A and Molecular Sieve 13X**  
(packed)



Molecular Sieve 5A 80/100 mesh  
1m x 1/8" x 2mm ID SilcoSmooth® tubing (cat.# 80440-800)



Molecular Sieve 13X 80/100 mesh  
2m x 1/8" x 2mm (ID) SilcoSmooth® tubing (cat.# 80439-800)