



USP
Pharmaceutical
Standards
See pages 527-528.

Organic Volatile Impurities (OVI) Analysis

Rtx®-G27 Column (fused silica with 5-meter Integra-Guard® guard column) (Crossbond® 5% diphenyl/95% dimethyl polysiloxane)

- Application-specific columns for residual solvents in pharmaceutical products. Meet all requirements of USP <467>.
- Analytical column with Integra-Guard® guard column eliminates connecting problems and leaks.
- Rtx®-G27 stable to 290 °C.

Some USP <467> methods require the use of a guard column. Our Integra-Guard® integrated guard column system makes it easy to comply.

ID	df	temp. limits	30-Meter with 5-Meter, 0.53mm ID Integra-Guard Guard Column
0.53mm	5.00µm	-60 to 270/290°C	10279-126

Rtx®-G43 Column (fused silica with 5-meter Integra-Guard® guard column) (Crossbond® 6% cyanopropylphenyl/94% dimethyl polysiloxane)

- Application-specific columns for residual solvents in pharmaceutical products. Meet all requirements of USP <467>.
- Analytical column with Integra-Guard® guard column eliminates connecting problems and leaks.
- Rtx®-G43 stable to 240 °C.

Some USP <467> methods require the use of a guard column. Our Integra-Guard® integrated guard column system makes it easy to comply.

ID	df	temp. limits	30-Meter with 5-Meter, 0.53mm ID Integra-Guard Guard Column
0.53mm	3.00µm	-20 to 240°C	16085-126

free literature

A Technical Guide for Static Headspace Analysis Using GC

Download your free copy from www.restek.com
lit. cat.# 59895A



Diane Thompson, Customer Service

Restek Customer Service

In the U.S.

Call: 800-356-1688 (ext. 3) or 814-353-1300 (ext. 3)

Monday–Friday 8:00 a.m.–6:00 p.m. ET

Fax: 814-353-1309—24-hours a day

Online: www.restek.com—24-hours a day

Outside the U.S.

Contact your Restek representative:
Refer to our list on pages 4-5 or visit our website at www.restek.com