



restek **innovation!**

Chlorinated Pesticides Analysis

Stx®-CLPesticides/Stx®-CLPesticides2

(proprietary Crossbond® phases)

- Application-specific columns for organochlorine pesticides and herbicides.
- Baseline separations in less than 10 minutes.
- Siltek® surface deactivation enhances responses for endrin, DDT, methoxychlor.
- Stable to 330 °C.

Many laboratories analyzing organochlorine pesticides struggle with breakdown and adsorption of endrin, DDT, and methoxychlor caused by active sites throughout the analytical system. Siltek® passivation technology enables these columns to offer unsurpassed inertness and the highest responses for active pesticides.

it's a **fact**

These columns are treated with Siltek® deactivation, which provides better responses for endrin, DDT, and methoxychlor.

Stx®-CLPesticides Columns (fused silica with Siltek® deactivation)

ID	df	temp. limits	15-Meter	30-Meter
0.25mm	0.25µm	-60 to 310/330°C	11540	11543
0.32mm	0.32µm	-60 to 310/330°C		11546
	0.50µm	-60 to 310/330°C	11541	11544

Stx®-CLPesticides2 Columns (fused silica with Siltek® deactivation)

ID	df	temp. limits	15-Meter	30-Meter
0.25mm	0.20µm	-60 to 310/330°C	11440	11443
0.32mm	0.25µm	-60 to 310/330°C	11441	11444

Stx®-CLPesticides Kits

(Note: Columns are not preconnected in these kits.)



Stx-CLPesticides Kit (0.25mm ID) cat.# 11190 (kit), \$1050 **SAVE \$100**

Includes:	cat.#	price
30m, 0.25mm ID, 0.25µm Stx-CLPesticides Column	11543	
30m, 0.25mm ID, 0.20µm Stx-CLPesticides2 Column	11443	
Universal Angled "Y" Press-Tight Connector, Deactivated	20403-261	
5m, 0.25mm ID Siltek Guard Column	10026	



Stx-CLPesticides Kit (0.32mm ID) cat.# 11193 (kit), \$1080 **SAVE \$135**

Includes:	cat.#	price
30m, 0.32mm ID, 0.32µm Stx-CLPesticides Column	11546	
30m, 0.32mm ID, 0.25µm Stx-CLPesticides2 Column	11444	
Universal Angled "Y" Press-Tight Connector, Deactivated	20403-261	
5m, 0.32mm ID Siltek Guard Column	10027	

Total cost if purchased separately \$1215

ordering **note**

Kits include Siltek® deactivated guard column.



Solutions For Your Environmental Analyses

Improved best-in-class GC columns • Standards • Industry experts at your service.

Visit us at www.restek.com/enviro