

LAMP GUIDE

ABI (APPLIED BIOSYSTEMS)

ANALYST	LH-13	D2	
ANALYST	LH-13X2	D2	
ANALYST	L-850	W	
BIOCAD 20, 60, 250, 700E, SPRINT, VISION	LH-40	D2	
INTEGRAL 700E, 250	LH-40	D2	
120A, 130A	LH-1/LH-4	D2	
120A, 130A	W		INQUIRE
270A-HT, 1000S	LH-1/LH-4	D2	
FS950	HG		INQUIRE
FS970, FS980	LH-1/LH-4	D2	
FS980	W		INQUIRE
PROSICE HT, CLC, 492, 494	LH-4	D2	
785	LH-4	D2	
SF757, 773, 783, 783A, 785A	LH-1/LH-4	D2	
SF757, 759A 773, 783, 783A, 785A	W		INQUIRE
SF759A	LH-4	D2	
SF769, SF770	LH-1	D2	

APPLIED BIOSYSTEMS /SEE ABI

ATI /SEE THERMO SCIENTIFIC

AVIV/PROTERION/RHEOMETRIC SCIENTIFIC/WYATT TECH CORP

14D			INQUIRE
62DS, 202, 215	L-995	XE	

BAS

FL45	LH-910	P/XE	
UV116A	LH-40	D2	
UV200	LH-40	D2	

BAUSCH & LOMB (B & L) / SEE THERMO SCIENTIFIC

AGILENT/HEWLETT-PACKARD

1037	L-610	W	
1040A, 1040M	LH-10	D2	
1046	L-890	P/XE	
1047A	L-610	W	
1084B	D2		INQUIRE
1050 SERIES			
79853A	L-295	D2	
79853C, 79854A	LH-10	D2	
G1306A	LH-10	D2	
79875A			INQUIRE
1090 SERIES			
79880A, 79881A	LH-10	D2	
SERIES I & II	LH-10	D2	
8450A	LH-16	D2	
8450A	LH-855	W	
8451A	LH-11	D2	
8452A	LH-12	D2	
8453A, 8453E			
8453A, 8453E	LH-13X2	D2	
8453A	L-850	W	
1100/1200 SERIES			
G1314A, G1314C	LH-14	D2	
G1314C	L-850	W	
G1315A,B & C; G1365A,B & C	LH-13	D2	
G1315A,B & C; G1365A,B & C	LH-13X2	D2	
G1315A,B & C; G1365A & B	L-850	W	
G1321A	L-870	P/XE	
G1362A		W	INQUIRE
HP3D/CE	LH-14	D2	

ALLTECH/GRACE

200, 201, 205	LH-40	D2	
200	LH-760	W	
305	LH-910	P/XE	

AMERSHAM

SEE GE HEALTHCARE/BIOSCIENCES/PHARMACIA/LKB

AMINCO /SEE THERMO SCIENTIFIC

AB-2		XE	INQUIRE
SLM 2000	L-280	D2	
SLM 3000	L-270	D2	
SLM 4800, 8000, 8100	L-855	XE	
SPF-500		XE	

BECKMAN/ALTEX

20, 24, 25, 26, 34, 35, 36	L-250	D2	
20, 24, 25, 26, 34, 35, 36	L-860	W	
25	L-200	D2	
153	L-885	HG	
160	L-710	HG	
160	L-715	ZN	
160	L-720	CD	
163	L-205	D2	
166	LHA-70	D2 Lamp Ass'y	
166	LH-70	Cartridge Only	
166		W	INQUIRE
168	LHA-75	D2 Lamp Ass'y	
168	LH-75	Cartridge Only	
5200, ACTA	L-250	D2	
5200, ACTA	L-865	W	
AD340	L-275	W	
	L-305	W	
DB, DBG, DBGT, DK, DU		D2	
DB, DBG, DK-2A	L-860	W	
DBGT		W	INQUIRE
DU2	L-60	W	
DU-6, 7, 8, 20, 30, 37, 40, 50	L-250	D2	
DU-6, 7, 8, 20, 30, 40, 50	L-840	W	
DU-60, 62, 63, 64, 65, 68, 70	L-835	W	
DU-60, 62, 63, 64, 65, 68	LH-71	D2	
DU-70	LH-73	D2	
DU-500, 520, 530	LH-110	D2	
DU-500, 520, 530	L-730	W	
DU-600, 640, 650, 800	LH-72	D2	
DU-600, 640, 650, 800	L-835	W	
DU-7000, 7400, 7500	LH-72	D2	
DU-7000, 7400, 7500	L-835	W	
P/ACE 2000, 2050	LH-77	D2	
P/ACE 2100, 2200, 2210	LH-77	D2	
P/ACE 5000/477170, 5010/477171	LH-76	D2	
P/ACE 5500/477172, 5510/477173	LH-76	D2	
P/ACE MDQ	LH-76	D2	
P/ACE 800	LH-76	D2	

BEIJING PURKINJE GENERAL INSTRUMENT CO./GENTECH

UV 1100, 9100, 1200, 1300, 1400	L-265	D2	
VIS-7220, 9100, 1200, 1300, 1400	L-165	W	



BIO-RAD			
760-1300	LH-25	D2	
760-1300	LHA-23	W	Lamp Assy
1305, 1306	L-290	D2	
1706, 1806A	L-295	D2	
1740		HG	INQUIRE
1790, 1801, BIO-DIMENSION	LH-40	D2	
1806B			INQUIRE
BIOFOCUS	LH-41	D2	
BIOLOGIC DUAL FLOW		HG	INQUIRE
EM-1 ECONO UV MONITOR		HG	
MICROSAMPLER 100	LH-41	D2	
REMEDI	LH-115	D2	
SMARTSPEC 3000	L-60	P/XE	
VARIANT		W	INQUIRE
VARIANT II		W	INQUIRE

BIOTAGE/DYAX			
HORIZON	L-735	HG	\$4
PARALLEX FLEX			INQUIRE
KILOPREP(100G)	LH-40	D2	
SP1	LH-60	D2	
FLASHMASTER II	LH-22	D2	

BIOTEK/ SEE KONTRON			

BISCHOFF			
LAMBDA 1000, 1010			INQUIRE
SP6, VISI-6			INQUIRE

BRANDTECH SCIENTIFIC			
WS2000			INQUIRE
WUV1101		XE	INQUIRE

BUCHI			
28062, 28063, 28067	LH-25	D2	

CAMAG			
TLC II SCANNER	L-205	D2	

CARLO ERBA/THERMO SCIENTIFIC			
	L-280	D2	

CARY/VARIAN			
1, 3, 4, 5, 100, 300, 400, 500	L-230	D2	
1, 3, 100, 300, 500	L-930	W	
400	L-110	W	
50		P/XE	INQUIRE
219	L-205	D2	
219	L-925	W	

CECIL			
CE 4300 WAVEQUEST			INQUIRE
CE 4200, 4201			INQUIRE
CE1200, 1220	LH-20	D2	
CE1020, 2020, 2040, 2041	L-210	D2	
CE 1020S			INQUIRE
CE 4500		XE	INQUIRE
CE 3000			INQUIRE
CE 1070			INQUIRE
CE 4000 - 9000 SERIES	L-210	D2	

COLE-PALMER			
01953-00	LH- 0	D2	

DIONEX/LC PACKING/GYNKOTEK			
39546, 39557			INQUIRE
AD20, PDA20		LH-5	D2
AD25, PDA100,PDA2996		LH-3	D2
PDA20			W
PDA100			W
CES I			INQUIRE
DSA-1		LH-40	D2
FDM-2			P/XE
RF1002/2000		L-800	XE
UDM-170S, 340S			INQUIRE
UDM-37182			INQUIRE
UVD 160, 170S, 170U, 320, 340S, 340U		LH-91	D2
VDM-1		L-200	D2
VDM-2		L-235	D2
UVD 3000			D2
UVD 3100, 3400			D2
UVD 3100, 3400			W

D-STAR			
DDW-10, DFL-10, DLC-20, DVW-10		LH-94	D2

DUPONT			
850		L-300	D2
8800		L-280	D2

DYAX /SEE BIOTAGE			

ELDEX			
MICROPRO			INQUIRE

EM SCIENCE /SEE HITACHI			

ERMA/ERC			
7511, 7512, 7515, 7522, 7530, 7535		L-610	W

ESA			
520, 522		LH-40	D2

FISHER/THERMO SCIENTIFIC			
20-0344		LH-40	D2

GBC			
AA		L-200	D2
CINTRA 5			INQUIRE
CINTRA 6, 10E, 20, 40, 916, 918, 920			INQUIRE
CINTRA 101,202,303,404,904			INQUIRE
CINTRA 906, 908, 911		L-200	D2
CINTRA 916		L-740	W
CINTRA 1155, 1255		LH-910	P/XE
CINTRA 1200			P/XE
CINTRA 1205, 1206, 1210		LH-40	D2
CINTRA 5100		LH-82	D2

GE HEALTHCARE/AMERSHAM/BIOSCIENCES/PHARMACIA/LKB			
2138, 2140, 2141			INQUIRE
2151		LH-20	D2
4050, 4051, 4053, 4054		L-210	D2
050		L- 35	W



GE HEALTHCARE (CONTINUED)			
GENEQUANT, GENEQUANT PRO			INQUIRE
ULTROSPEC II, III, PLUS	L-210	D2	
ULTROSPEC II, III, PLUS	L-835	W	
ULTROSPEC 100 PRO		W	
ULTROSPEC 500	L-935	W	
NOVASPEC II	L-835	W	
NOVA SPEC III		W	
NOVA SPEC PLUS		W	
ULTROSPEC 1000E, 2000, 3000	LH-125	D2	
ULTROSPEC 3300, 4000, 4300	LH-125	D2	
ULTROSPEC 1100 PRO, 3000 PRO	LH-125	D2	
ULTROSPEC 1000E, 2000, 3000	L-935	W	
ULTROSPEC 3300, 4000, 4300	L-935	W	
ULTROSPEC 1100 PRO, 3000 PRO	L-935	W	
ULTROSPEC 2100 PRO	L-640	P/XE	
ULTROSPEC 2158			INQUIRE
ULTROSPEC 2158		W	INQUIRE
ULTROSPEC 3100	L-640	P/XE	
UV 900	L-990	P/XE	
UVIS 920		HG	INQUIRE
UPC900		HG	INQUIRE
AKTA DESIGN SYSTEMS			
BASIC, EXPLORER, PILOT, PURIFIER	L-990	P/XE	
FPLC		HG	INQUIRE
PROCESS			INQUIRE

GILFORD			
200, 220, 222, 230, 240, 250, 252, 260	L-205	D2	
2000, 2400, 2600	L-205	D2	
200, 222, 230, 240, 250, 252, 260	L-875	W	
2000, 2400, 2600	L-875	W	
RESPONSE I	L-210	D2	
RESPONSE II			
RESPONSE I	L-900	W	
RESPONSE II			
STASAR I, II, III, SYVA S-3	L-705	W	
FLUORO IV	L-801	XE	

GILSON			
111			INQUIRE
115, 116, 117	LH-150	D2	
118, 119, 151, 152, 155, 156	LH-151	D2	
122	L-800	XE	
160			INQUIRE
170			INQUIRE
170		W	INQUIRE
HOLOCHROME			INQUIRE
SPECTROCHROME		W	INQUIRE

GROTON			
20/20			INQUIRE
FD100	L-890	P/XE	

GYNKOTEK /SEE DIONEX			

HACH			
DR/4000U	LH-110	D2	\$4
DR/4000U, DR/4000V	L-730	W	
DR/5000			INQUIRE
DR/5000		W	INQUIRE
000, 100	L-730	W	

HEWLETT-PACKARD /SEE AGILENT			

HITACHI			
100-10, 20, 30, 40, 60, 80	L-280	D2	
100-10, 20, 30, 40, 60, 80	L-805	W	
101, 102, 111	L-430	D2	
101, 102	L-805	W	
124	L-210	D2	
124	L-805	W	
150-20	LH-51	D2	
150-20	L-805	W	
181			INQUIRE
181, 191	L-805	W	
200, 300, 330, 340	LH-51	D2	
655A-22, 23	LH-51	D2	
F-1000, 1050, 1080, 1200	L-800	XE	
F-2000, 2500, 3000, 3010	L-800	XE	
F-4000, 4010, 4500	L-800	XE	
F-4500			INQUIRE
F-7000	L-800	XE	
GENE SPEC I, III	LH-56		
GENE SPEC V			INQUIRE
L-2400, 2420, 2450	LH-56	D2	
L-2400, 2420, 2450, 2455	L-835	W	
L-2480, 2485	L- 00	XE	
L-2500, 3000, 3100, 3200	LH-51	D2	
L-4000	LH-56	D2	
L-4200	LH-51	D2	
L-4200	L-815	W	
L-4250	LH-51	D2	
	L-815	W	
L-4500	LH-56	D2	
L-7000, 7400, 7420, 7450A	LH-56	D2	
L-7300, 7455	LH-56	D2	
L-7420	L-815	W	
L-7480, 7485, 7490	L-800	XE	
SERIES 100	L-805	W	
U-0080D	LH-56	D2	
U-1100, 1500, 1800	LH-51	D2	
U-1100, 1800	L-815	W	
U-1500, 3210	LH-55	W	
U-1900	LH-56	D2	
U-1900	L-815	W	
U-2000, 2001, 2010, 2810	LH-51	D2	
U-2000, 2001, 2010	L-815	W	
U-2900, 2910	LH-56	D2	
U-2900, 2910	L-815	W	
U-3000, 3010, 3110	LH-51	D2	
U-3210, 3300, 3310	LH-51	D2	
U-3000, 3010, 3110	L-815	W	
U-3300, 3310	L-815	W	
U-3400, 3410, 3501	LH-51	D2	
U-4001, 4100	LH-51	D2	
U-3410, 3501, 4001	LH-55	W	
U-3500	LH-56	D2	
U-7000			INQUIRE

HNU			
103C	L-750	PID	\$

HORIBA JOBIN YVON /SEE JY			

IBM			
9420	L-200	D2	
9523, 9563, 9565	L-205	D2	



IONICS BUSINESS GROUP /SEE SIEVERS

ISCO/TELEDYNE			
205	LH-40	D2	
205		W	INQUIRE
3850	L-205		
COMBIFLASH BIO OPTIX 10	LH-61		
COMBIFLASH COMPANION, 4X, XL, RF	LH-61		
COMBIFLASH SYSTEM SQ16X	LH-61		
V4/1840	L-205	D2	

KNAUER/SONNTEK			
200 WC		LH-2	HG
2000, 2001, 2500, 2501 WC		LH-25	D2
2600, 2650, 2700 WC		LH-25	D2
2001, 2500, 2501 WC		LHA-23	W
2500 SL, 2600 SL		LH-25	D2
2600, 2700 SL		LHA-23	W
2800 SL		LH-26	D2
2800 SL			W INQUIRE
7100		LH-21	HG
8700, 9700		LH-20	D2
D871		LH-20	D2
D2540		LH-21	HG

JASCO			
870, 875	L-295	D2	
870		W	INQUIRE
7709, 7800			INQUIRE
7800		W	INQUIRE
FP-750, 777	L-800	XE	
FP-820, 821, 920	L-800	XE	
FP-1520, 1570	L-800	XE	
FP-1575, 2020	L-800	XE	
FP-3020			INQUIRE
FP-3120	L-800	XE	
FP-6200, 6300, 6500, 6600	L-800	XE	
MD-910, 915, 1510, 1515, 2010			INQUIRE
MD-910, 915, 1510, 1515, 2011		W	INQUIRE
MD-2015			INQUIRE
OR-2090		EX-HG	INQUIRE
CD-2095			INQUIRE
UV-970A, 975A	L-295	D2	
UV-970B, 975B, 975C			INQUIRE
UV-1570, 1575, 2070, 2075			INQUIRE
UV-2070		W	INQUIRE
UV-2077			INQUIRE
UVIDEC 100, 100-II THRU -VI	L-200	D2	
UVIDEC 320, 340, 420, 430	L-200	D2	
UVIDEC 505, 510, 610, 650, 710	L-200	D2	
V-530	L-265	D2	
V-550	L-835	W	
V-560, 570, 600, 630 BIO			INQUIRE
V-7000, 7200, 7300			INQUIRE

KONTRON/BIOTEK /NORTH STAR SCIENTIFIC			
430			INQUIRE
440		LHA-34	D2
720, 810, 820, 860		L-205	D2
820, 860		L-110	W
900, 930, 940		L-205	D2
900, 930, 941		L-110	W
UVIKON		L-205	D2
UVIKON XL, XS			INQUIRE
ALL SPECTROPHOTOMETERS		L-205	D2

KRATOS /SEE ABI/PERKIN-ELMER/SCHOEFFEL

LABALLIANCE/ SCIENTIFIC SYSTEMS (SSI)			
200, 201		LH-40	D2
500, 504, 506		LH-40	D2
V520, 525		LH-40	D2

LC PACKINGS/DIONEX			
ULT-UVD, UVD-3000		LH-25	D2
ULT-UVD		LHA-23	W
ULT-UVD		LH-23	W INQUIRE

LDC, LINEAR /SEE THERMO SCIENTIFIC

LKB /SEE GE HEALTHCARE			
4050, 4051, 4054		L-210	D2
4050		L-835	W

MCPHERSON			
FL-750/A & B		LH-101	XE
FL-750/A & B		LH-100	D2
FL-750/A & B		LH-102	XE-HG

MERCK / SEE HITACHI

MICHROM BIORESOURCES			
UV-200		LH-40	D2
MD1 THRU MD4		LH-25	D2

MILTON-ROY			
SEE THERMO SCIENTIFIC, BAUSCH&LOMB			

JENWAY			
6100		W	INQUIRE
6105, 6405, 6505	LH-650	D2	
6105, 6300, 6400, 6405, 6500, 6505	L-835	W	
6200, 6305		XE	INQUIRE
AQUANOVA		XE	INQUIRE
		XE	INQUIRE

JY (JOBIN YVON)/SPEX/HORIBA			
1934D		P/XE	INQUIRE
FLUOROLOG-2/F111AI, 112AI		XE	INQUIRE
FLUOROLOG-2/F111XI, 112XI	L-855	XE	
FLUOROLOG-2/F1T11I	L-855	XE	
FLUOROLOG-2/F212I, F222I	L-855	XE	
FLUOROLOG-3	L-855	XE	
FLUOROMAX-2, 3, 4	L-955	XE	
SYSTEM 5000	L-60	P/XE	



MOLECULAR DEVICES		
FLEX STATION		P/XE
SPECTRA MAX 340	L-65	P/XE
SPECTRA MAX 190, 340PC	L-60	P/XE
SPECTRA MAX M2/M2E, M5/M5E	L-870	P/XE
SPECTRA MAX ANALYST GT	L-760	P/XE
SPECTRA MAX GEMINI, XS, XPS, EM	L-760	P/XE
SPECTRA MAX 250	L-760	P/XE
SPECTRA MAX PLUS, 384 PLUS	L-60	P/XE
UV MAX, THERMO MAX		W
VMAX		W

NICOLET /SEE THERMO SCIENTIFIC

OLIS		
RSM-1000	L-1000	XE

PERKIN-ELMER/EG&G

HPLC UV Detectors		
785A	LH-98	D2
INTEGRAL 4000	L-210	D2
LC-55, 55B, 65, 75	L-210	D2
LC-85B, 95, 135, 235	L-210	D2
LC-90		INQUIRE
LC-200	LH-98	D2
LC-200		W INQUIRE
LC-240		P/XE INQUIRE
LC-250, 290, 295		INQUIRE
LC-480	LH-92	D2
SERIES 200 /DAD, PHASE I, PHASE II	LH-97	D2
SERIES 200 /DAD	L-725	W
SERIES 200 /FLUORO	L-800	D2
SERIES 200 /UV-VIS/785A	LH-98	D2
SERIES 200 /UV-VIS/785A Fluorometers		W INQUIRE
SERIES 203, 204, 512, 650	L-800	XE
LS-3, 4, 5		P/XE INQUIRE
LS-30, 40, 50		P/XE INQUIRE
LC-240		INQUIRE
MPF-2A, 3, 4	L-800	XE
MPF-43, 43A, 44, 44A, 44B	L-800	XE
UV-Visible Spectrophotometers		
46, 55, 57	L-200	D2
200		INQUIRE
220, 320, 330, 340	LH-51	D2
550, 551, 552, 555	L-240	D2
550S, 551S, 554	L-400	D2
552A, 553, 559, 559A	L-200	D2
556, 557	LH-51	D2
570, 571, 572, 575, 576	L-200	D2
3840	L-200	D2
BIO	LH-92	D2
BIO	LH- 5	W
EZ150, 201, 210	LH-51	D2
EZ150, 201, 210	L- 0	W
LAMBDA 1	L-300	D2
LAMBDA 1	LH-845	W
LAMBDA 2, 5, 7, 9, 10, 11, 12	LH-92	D2 \$
LAMBDA 2, 3B, 5, 7, 9, 10, 11, 12	LH- 5	W
LAMBDA 3, 3B	L-300	D2
LAMBDA , 6	LH-93	D2
LAMBDA , 6	LH- 5	W
LAMBDA 14 THRU 19	LH-92	D2
LAMBDA 14 THRU 19	LH- 5	W \$
LAMBDA 20, 25, 35, 0, 5	LH-92	D2

PERKIN-ELMER (CONTINUED)			
LAMBDA 20, 40	LH-845	W	
LAMBDA 650, 800, 850, 900, 950	LH-92	W	
LAMBDA 650, 800, 850, 900, 950	LH-845	W	
LAMBDA XLS			INQUIRE
Atomic Absorption Spectrometers			
100, 300 ANALYSIS	L-335	D2	
200, 400 ANALYSIS	L-410	D2	
272, 280	L-310	D2	
303, 305, 305A, 305B, 306			INQUIRE
360, 370, 372, 373, 380	L-310	D2	
402, 403			INQUIRE
410, 420, 430	L-410	D2	
460, 560	L-310	D2	
480, 503, 555, 575			INQUIRE
560	L-310	D2	
600, 700, 800 ANALYSIS	L-405	D2	
603, 605, 703			INQUIRE
1100, 1100B, 2100	L- 05	D2	
2280, 2380	L-310	D2	
2380	L-310	D2	
3030B	L-310		INQUIRE
3100, 3110, 3300 ANALYSIS	L-335	D2	
4000 ANALYSIS	L-310	D2	
4100	L-405	D2	
5000, 5100 ANALYSIS			
Polarimeters			
241, 241 MC	L-600	NA	
243, 341, 343	L-600	NA	
Plate Readers			
2100, 2103 ENVISION		P/XE	INQUIRE

PERSEPTIVE BIOSYSTEMS/SEE ABI

PHARMACIA /SEE GE HEALTHCARE

PHOTON TECHNOLOGY INT'L (PTI)

QM-2/2005, 4/2005, 6/2005, 7/2005, 8/2005	L-970	XE
QM-14/2005	L-970	XE

POSTNOVA

S3210	LH-13	D2
S3240	LH-13	D2

ROSEMOUNT ANALYTICAL

90	L-330	D2
----	-------	----

SCIENTIFIC SYSTEMS /SEE LABALLIANCE

SEDEX

55, 65, 75, 75C	L-935	W
55LT, 60LT, 75LT		
75		

SEQUOIA TURNER

340, 360, 690		W	\$
---------------	--	---	----



SHIMADZU			
HPLC UV Detectors			
LC2010AHT, LC2010CHT	LH-	D2	25.00
SP-4	L-305	D2	\$295.00
SPD-2A, 2AV, 6A, 6AV	L-305	D2	\$295.00
SPD-6AV	L-100	W	\$35.00
SPD-7A, 7AV	L-305	D2	\$295.00
SPD-10A, 10AV, 10AVP, 10AVVP	LH-82	D2	\$450.00
SPD-10AV, 10AVVP	L- 30	W	15 .00
SPD20A, 20AV	LH- 2	D2	50.00
SPD20AV	L-830	W	\$158.00
SPD-M1A			INQUIRE
SPD-M6A, 10A, 10AV	LH-81	D2	\$450.00
SPD-M6A, 10A	L-935	W	\$28.00
SPD-M10AV, -M10AVP	L-245	W	\$158.00
SPD-M10AVP	LH-88	D2	\$450.00
SPD-M20AV	LH-87	D2	
SPD-M20AV	L-830	W	
Fluorometers			
RF-10A, 10AXL	L-800	XE	
RF-10AVP			INQUIRE
RF-530	L-975	XE	
RF-535, 540, 551	L-800	XE	
RF-551	L-801	XE	
RF-1501, 5000, 5301	L-800	XE	
UV-Visible Spectrophotometers			
BIOSPEC 1601	L-265	D2	
BIOSPEC 1601	L-110	W	
CS-920, 930	L-265	D2	
CS-930	L-100	W	
CS-9000U, 9301	L-975	XE	
CS-9301	L-265	D2	
MULTISPEC 1501	L-265	D2	
MULTISPEC 1501	L-935	W	
SERIES 160	L-110	W	
UV-120, 120-02, 140-02	L-265	D2	
UV-120, 160U	L-110	W	
UV-150, 150-02, 160U, 163	L-265	D2	
UV-240, 250, 260, 265	L-265	D2	
UV-240, 260, 265	L-110	W	
UV-1201, 1201PC, 1240	L-265	D2	
UV-1201, 1240, 1240V	L-245	W	
UV-1501, 1600, 1601, 1601PC, 1603, 1650, 1700, 1800	L-265	D2	
UV-1501, 1600, 1601, 1601PC, 1603, 1650, 1700, 1800	L-110	W	
UV-2100, 2101, 2401PC, 2501PC	L-265	D2	
UV-2501PC, 3101PC	L-110	W	
UV-3000			INQUIRE
UV-3100, 3100PC, 3101PC	L-265	D2	
UV-Visible Spectrophotometers			
UV2101PC, UV-2401, 2450, 2550	L-110	W	
UV3600	L-110	W	
Atomic Absorption Spectrometers			
AA-625, 630, 640, 660	L-265	D2	
AA-670, 680U	L-265	D2	
AA6200	L-265	D2	
AA-6501, 6601, 6701, 6800			INQUIRE

SHODEX			
RI71, 72, 73, 74			
RI101			

SIEVERS/IONICS INSTR. BUSINESS GROUP/GE ANALYTICAL INSTR.			
400, 800	L-700	HG	
900			INQUIRE
PPT	L-605	HG	

SPECTRA-PHYSICS /SEE THERMO SCIENTIFIC			
SPECTRONICS/ SEE THERMO SCIENTIFIC			
SPEX /SEE JY (JOBIN YVON)			
TEKMAR-DOHRMAN /SEE TELEDYNE			
TELEDYNE / TEKMAR INSTRUMENTS			
PHOENIX 8000			INQUIRE
TOX ROUND GRID LAMP			INQUIRE
OZONE GENERATOR LAMP 7			INQUIRE
TECAN			
SPECTRO FLUR +		P/XE	\$595.00
GENIOS		P/XE	\$595.00

THERMO SCIENTIFIC			
ATI/UNICAM/NICOLET			
ATI/240, 250	LHA-34	D2	
ATI/240, 250	LH-34	D2	
HELIOS ALPHA, BETA, DELTA, GAMMA	LH-54	D2	
HELIOS ALPHA, BETA, DELTA, GAMMA	L-910	W	
NICOLET EVOLUTION 100	LH-54	D2	
NICOLET EVOLUTION 100	L-910	W	
NICOLET EVOLUTION 300		P/XE	INQUIRE
NICOLET EVOLUTION 500, 600	LH-54	D2	
NICOLET EVOLUTION 500, 600	L-940	W	
HELIOS EPSILON	L-915	W	
UNICAM UV 1-100, 2-100, 2-200	LH-54	D2	
UNICAM UV 1-100	L-910	W	
UNICAM UV 2-100, 2-200	L-940	W	
UNICAM UV 2-300, 2-400, 3-100	LH-54	D2	
UNICAM UV 2-300, 2-400, 3-100	L-940	W	
UNICAM UV 3-200, 4-100, 4-200	LH-54	D2	
UNICAM UV 3-200, 4-100, 4-200	L-940	W	
UNICAM UV 4-500	LH-54	D2	
UNICAM UV 4-500	L-940	W	
UNICAM UV 300, 500, 550	LH-54	D2	
UNICAM UV 300, 500, 550	L-940	W	
UNICAM SP6550	L-245	W	
UNICAM SP8800	L-255	D2	
BAUSCH & LOMB (B & L)			
AQUAMATE, AQUAMATE UV	L-910	W	
AQUAMATE UV	LH-54	D2	
BIOMATE 3		P/XE	INQUIRE
BIOMATE 5	L-910	W	
DE 50	L-225	D2	
GENESYS 2, 5	LH-105	D2	
GENESYS 2, 5	L-915	W	
GENESYS 6		P/XE	INQUIRE
GENESYS 8	LH-54	D2	
GENESYS 8	L-910	W	
GENESYS 10, 20	L-915	W	
GENESYS 10 UV		P/XE	INQUIRE
SERIES II		XE	INQUIRE
SPECTRONIC 20	L-905	W	
SPECTRONIC 21	L-945	W	
SPECTRONIC 21 DUV	L-225	D2	
EPSILON	L-915	W	
SPECTRONIC 70, 100, 700, 701, 710	L-920	W	
SPECTRONIC 120			INQUIRE
SPECTRONIC 200, 601, 700, 701, 710	L-225	D2	



THERMO SCIENTIFIC (CONTINUED)			
SPECTRONIC 301, 401, 501, 601	L-615 *	W	
SPECTRONIC 301, 401, 501, 601	L-915 **	W	
SPECTRONIC 501, 601			INQUIRE
SPECTRONIC 1001, 1201, 2000	L-225	D2	
SPECTRONIC 1001, 1201	L-810	W	
SPECTRONIC 2000	L-920	W	
SPECTRONIC 3000	L-270	D2	
SPECTRONIC 3000		W	INQUIRE
ENVIRONMENTAL INSTRUMENTS			
43A, CT, CH	L-40	P/XE	
LDC			
120, 1204A			INQUIRE
SC100-234	LH-40	D2	
SC100-234	LH-760	W	
SM I, II, III, D	L-205	D2	
SM 3000, 3100, 4000	L-205	D2	
SM 3100X, 4100			INQUIRE
SM 3200			INQUIRE
SM 5000	LH-47	D2	
UV100-411, UV150-511	LH-40	D2	
UV100-411, UV150-511	LH-760	W	
LINEAR			
100			INQUIRE
200, 203, 204, 205, 206PHD	LH-40	D2	
200, 201, 203, 204, 205	LH-760	W	
304, 305	LH-910	P/XE	
SPECTRA-PHYSICS/FINNIGAN			
CE500, 1000, 2000	LH-40	D2	
CE500, 1000, 2000	LH-760	W	
FL2000, 3000	LH-910	P/XE	
SC100, 200, SF100, 200	LH-40	D2	
SC100, 200, SF100, 200	LH-760	W	
SF770, 970	LHA-1	D2	
SF770, 970	LH-1	D2	
SP8400, 8440, 8440XR	L-215	D2	
SP8450, SP8490	L-215	D2	
SP8480XR	LHA-3	D2	
SP8480XR	LH-1	D2	
SP8773XR	LHA-2	D2	
SP8773XR	LH-1	D2	
UV100, 150, 1000, 2000, 3000	LH-40	D2	
UV100, 150, 1000, 2000, 3000	LH-760	W	
UV3000HR	LH-40	D2	
UV6000LP	LH-46	D2	
UV6000LP	L-765	W	
TURNER/BARNSTEAD			
110, 111, 112	L-185	HG	
SP-830+, 850, 870+		W	
SP-890/SM110245			INQUIRE
* L-615 ROUNDED TOP, ** L-915 POINTED TOP			

TRANSGENOMIC			
FDX-99-2100	L-800	XE	
WAVE 3500, 3500A, 3500HT	LH-56	D2	\$395.00
WAVE HS	L-800	XE	
WAVE MD		HG	INQUIRE
WAVE OLIGO	LH-56	D2	

TURNER/BARNSTEAD /SEE THERMO SCIENTIFIC

UNICAM /SEE THERMO SCIENTIFIC

VARIAN/CARY/RAININ			
1, 3, 4, 5, 100, 300, 400, 500, 2300	L-230	D2	
1, 3, 100, 00, 500	L-930	W	
50		P/XE	INQUIRE
7/ 700	LH-20	D2	
14, 15, 16, 17, 210, 219, 220	L-230	D2	
14, 15, 16, 17, 210, 219, 220	L-930	W	
300	L-110	W	
2050	L-290	D2	
2070, 9050, 9060, 9065			INQUIRE
2550	L-295	D2	
8700	LH-20	D2	
9070, 9075	LH-910	P/XE	
AA200, 275, 300, 400	L-230	D2	
AA600, 640, 775, 800, 875	L-230	D2	
AA975, 1275, 1475	L-230	D2	
AA3000	L-230	D2	
DMS90			INQUIRE
ECLIPSE			INQUIRE
FL-1, 2	LH-910	P/XE	
LC-5000, 5500	L-305	D2	
PDA-1	LH-81	D2	
PDA-2	LH-82	D2	
PDA-2		W	INQUIRE
PRO STAR, PRO STAR 320	LH-950	P/XE	
PRO STAR 310	L-305	D2	
PRO STAR 325, 335	LH-325	D2	
PRO STAR 330	LH-56	D2	
PRO STAR 340, 345	LH-40	D2	
PRO STAR 340, 345	LH-760	W	
PRO STAR 360	LH-910	P/XE	
PRO STAR 363	L800	XE	
PRO STAR 364		P/XE	INQUIRE
SPECTRA A10, 20, 30, 40	L-230	D2	
STAR 9050	L-305	D2	
UV-1, 1A	LH-950	P/XE	
UV50, 100, 200	L-305	D2	
UV-C, DII, M	LH-40	D2	
VISTA 5500	L-305	D2	

WATERS			
401, 403, 404			ADVISE
420	L-770	HG	
440, 441	L-820	HG	
441	L-720	CD	
441	L-715	ZN	
450	LHA-1	D2 Lamp Ass'y	
450	LH-1	D2 Replacement	
470, 474	L-800	XE	
480, 481	L-205	D2	
484	LH-30	D2	
486	LH-31	D2	
490E	L-890	P/XE	
990, 991, 994	LHA-34	D2 Lamp Ass'y	
990, 991, 994	LH-34	D2 Replacement	
2996, 996	LH-33	D2	
2475	LH-925	XE	
2487, 2488	LH-36X2	D2	
ACQUITY UPLC FLR DETECTOR		HG/XE	ADVISE
ACQUITY UPLC PDA DETECTOR		PAPF	
ACQUITY UPLC TUV DETECTOR		PAPF	
LC MODULE 1	LH-31	D2	

ZEISS
MCS 501 **INQUIRE**



