



MULTI VORTEX-GENIE®

Large Capacity Multi Vessel Vortex Mixer



Vortex Anything and Everything that only a Vortex-Genie® can - Only More!

Another revolutionary product from the makers of the Vortex-Genie!

- ⦿ Same true vortexing as with the original Vortex-Genie® 2.
- ⦿ Holds multiple test tubes, vials, or cylinders of various sizes (racks included) or up to six 250mL vessels, (more for smaller vessels) or a combination of tubes and vessels such as Erlenmeyer flasks.
- ⦿ Designed and built "Genie™" strong, rugged and reliable.
- ⦿ Digital control and display of speed and time for reproducibility and accuracy, with true speed display.
- ⦿ Gentle ramp up to speed, so open containers will not spill.
- ⦿ Use for gentle, orbital-style shaking or more vigorous vortexing.
- ⦿ Durable and heavy zinc die casting with specially designed feet prevents "walking".
- ⦿ No additional holders or accessories necessary. Racks are included.
- ⦿ May be used in cold rooms or incubators.

Ordering Information

Catalog No.	Description
SI-M236	Multi Vortex-Genie, 120V, 60Hz
SI-M246	Multi Vortex-Genie, 230V, 50Hz, No Plug
SI-M256	Multi Vortex-Genie, 230V, 50Hz, Euro Plug
SI-M266	Multi Vortex-Genie, 230V, 50Hz, British Plug
SI-M276	Multi Vortex-Genie, 230V, 50Hz, Swiss Plug
SI-M296	Multi Vortex-Genie, 230V, 50Hz, Aust. Plug
SI-M286	Multi Vortex-Genie, 100V, 50/60Hz

Specifications

Supplied (Capacity)	Two 12-Pack Multi-Tube Racks holding up to 24 tubes of any material up to 29mm dia., or up to six 250mL vessels.
Speed	250 - 1000RPM (up to 950RPM for 50Hz models)
Vortexing Orbit	4mm dia. (same as Vortex-Genie 2)
Timer	1-99 minutes or continuous
Dimensions (D x W x H)	310 x 270 x 195mm (12.2 x 10.6 x 7.7in)
Weight	11.3 kgs (25 lbs)



Home of the Multi-Task Genie™ Family
 70 Orville Drive, Bohemia, NY 11716 U.S.A.
 Phone: (631) 567-4700 • Fax: (631) 567-5896
 Toll Free: 888-850-6208
 Customer Service: custsvc@scientificindustries.com
www.scientificindustries.com